ivinito^{iq}

ONE-BODY,
PRESS-OVER-METAL
PRESS-OVER-ZIRCON

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

GC Initial IQ • One-Body, Press-over-Metal Opaque Layering

GC Initial IQ • One-Body, Press-over-Metal Opaquerschichtung

Cast or CAD/CAM made metal frameworks are finished using tungsten carbide metal burs or ceramic-bonded grinding instruments and oxidized according to alloy manufacturers instructions.

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Apply the first Opaque layer (wash) with the normal GC Initial MC Paste Opaques in a thin layer using a flat brush. Before firing fluorescent crystals are applied on the Opaque surface

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After firing these Fluo Crystals are the base for a very important equal layer of the 2nd Opaque firing.

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Fire 2nd Opaque layer as usual, at critical areas the coloured Opaque Modifiers should be used. The Opaque layer should have a slightly shiny surface. Beside the colours A1, B1, C1, a one shade darker Opaker than the requested colour is being applied.

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GC Initial IQ • One-Body, Press-over-Zircon Power Frame Modifiers

The pre-coloured Zirconium dioxide framework.

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The pre-coloured Zirconium frame structures can be changed in colour and fluorescence by applying a thin wash layer of Power Frame Modifiers.

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It is unfortunately impossible to offer a general colour combination chart for Power Frame Modifiers with the 16 Vitapan Classic® shades as the existing Zr frame colouring liquids do vary too much from manufacturer to manufacturer.

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*GC Initial IQ • Power Frame Modifiers in the posterior regions.
Aesthetical critical area's should be supported by extra colour from the Power Frame Modifiers.*

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Full contour waxing up Voll anatomische Wachsmodellation

After the Opaque firings, the restoration is completely waxed up (anatomical & functional).

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With Press-over-Metal, a ceramic shoulder is aesthetically inevitable. With Press-over-Zircon, a ceramic shoulder is not recommendable. The outline area of a Zircon based restoration should always be supported by the Zircon frame structure.

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Full wax up with the aid of a pre-fab wax-veneer system.

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Minimum thickness, labial/buccal for pressed objects is 0,8mm to reach an acceptable aesthetical endresult.

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Full contour waxing up Voll anatomische Wachsmodellation

For excellent full anatomic, diagnostic wax up, it is recommendable to use one, on the market available, pre-fab wax veneer system. The system used in this Technical Manual is “Anteriores” from MDT. Wichnalek. This system combines resin models with the corresponding wax veneers and the clinical “Anteriores - Picture Gallery” book from Dr. med. Dent. Jan Hajtó (München – Germany).

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www.teamwork-media.de

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Sprueing of the modellation:

The wax sprues (3.0-3.5mmØ) should be app. 5-6 mm in length and should not taper towards the pattern (Rounded out – no sharp edges, consider the direction of flow of the ceramic material: same flow direction).

The patterns are waxed carefully onto the ring base former in an angle of 45° (Distance between the individual objects: min. 3mm).

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
Remark: Multi-unit bridge constructions should be made with a crossbar instead of direct sprueing.

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Wax weight: Please weight the objects before investing. To get the right wax-weight, the frame structure has to be weighed before and after the anatomical wax modellation.

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The investing is carried out with GC MultiPressVest, specially developed, fine grained phosphate bonded investment (For details please refer to the special instruction for use).

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Mix GC MultiPressVest acc. the instructions for use:
premix powder/liquid by hand with a spatula until a uniform mixture.

Mixing ratio:

Small ring size: 100g powder – 22ml liquid

Medium ring size: 200g powder – 44ml liquid

Large ring size: 300g powder – 66ml liquid

Mix for 60 sec. under vacuum (320 – 420 rpm)

After filling the investment ring carefully, put on the mould base former turning it gently (Note: the investment needs to extrude easily through the opening). After setting remove the mould base by turning carefully. Push the investment cylinder out of the silicone ring. Smooth the underside by using a plaster knife (90° angle, stable position).

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Notes:

- Check the correct fit of the silicone cylinder before investing.
- The base of the cylinder must be absolutely flat (90° angle), thus the mould stands perfectly in an upright, stable position in the injection furnace.
- Investment residues may not enter the sprues. Clean them carefully before preheating.
- Investment materials contain quartz powder. Avoid the inhalation of dust
- Please observe the corresponding Instruction for Use for GC MultiPressVest.

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Preheating:

A.: Quick heating method:

After 20 min. put the investment (the funnel pointing downwards) centrally in the furnace, pre-heated to 850°C (Pre-heating: 60min).

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If there is more than 1 ring in the furnace, each holding time should be extended by 10 min per add. cylinder.

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- GC Initial IQ Press Pellets should not be preheated.
- Alox-press plungers can be put straight at 850°C in the oven.
- One way press plungers do not need to be pre-heated.
- All Alox-press plungers need to be clean before pressing (Use glassbeads for sandblasting).
- Attention: short timing and distance between preheating furnace and pressing is needed!

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Pressing programs for various calibrated press furnaces

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All data are general references of in house tests of calibrated press furnaces. Please calibrate your press furnaces regularly and test the final temperature. In case of too low pressing temperatures, the CTE-value of press ceramic may change, resulting in cracking and chipping.

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After cooling, please mark the length of the plunger on the investment ring and cut separate along the length of the mark using a suitable cut-off disc (sintered diamond disc).

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Afterwards carefully break apart the 2 sections.

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Roughly sandblast the injected patterns with glass beads (4 bar, 50 μ m), after this clean the objects carefully and completely (2 bar, 50 μ m).

Attention: Do not use Aluminium oxide!

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Devesting and finishing

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Cut of the sprues carefully using a suitable cut-off disc (fine diamond disc, no pressure and low rotation). The ceramic objects are prepared according to the technique used (fine diamonds, low RPM, avoid heating up to prevent cracking). The contact areas metal frame / ceramic (e.g. ceramic shoulder) and the Zircon frame structure / ceramic have to be treated with caution. In case of small failures, a correction powder of each ingot colour is available in the sets.

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Initial IQ • Lustre Pastes Application

Initial IQ • Lasurenpastenanwendung

Before applying the Lustre Coating, the ceramic restoration is shaped and contoured using adjusted diamond burs. The finished pressed-over Crowns & Bridges are sandblasted with 50 microns Alu-Oxid at 1,5 bar pressure.

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The Lustre Paste Neutral is coated on the entire surface of the sandblasted restoration. This coating is clearly thicker than what we know from a normal Glaze firing.

Note: these Lustre Pastes may not come in contact with water.

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The requested Vitapan Classic® Shade is created by a more or less intensive application of the Lustre Pastes "Body A – B". This process, which only takes 2 minutes, is everything which is necessary to create Vitapan Classic® shades.

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Initial IQ • Lustre Pastes Application

Initial IQ • Lasurenpastenanwendung

In order to achieve the higher grey value of the incisal area with darker teeth (A3.5; A4; C3; C4), 2 grey Effect Shades can be used: L-3 (Light Grey) and L-4 (Dark Grey).

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For eventual further individualisation, you can use the 2 whitish Lustre Pastes: L-1 (Vanilla) and L-2 (White) for brightening the restorations.

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2 bluish Effect Lustre Pastes: L-5 (Light Blue) and L-6 (Dark Blue) for a bluish opalescent appearance in the incisal area.

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In order to achieve an even higher „Transparency“ the 2 greyish Effect Lustre Pastes L-3/L-4 can be used as well.

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All GC Initial MC INvivo Stains can be used to achieve even more possibilities of individualisation.

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Initial IQ • Lustre Pastes Application

Initial IQ • Lasurenpastenanwendung

The requested surface texture / smoothness of the Lustre Paste layer can be obtained by soft vibrating or condensing the restoration.

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
Final result after 1 firing. The procedure for metal supported & ZrO₂ supported restorations is identical. In case of Press-over-Zircon restorations, the Lustre Pastes layering should be fired in 2 thin layers instead of 1.

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

Initial IQ Initial IQ

A short overview of the capacity of the GC Initial IQ • One-Body, Press-over-Systems.

All Press-over-Metal restorations only had one Lustre Paste firing, the Press-over-Zircon restorations had 2 firings.

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All technical / physical values provided here, refer to in-house testing following EN ISO 9693 : 2000.

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GC Initial IQ • One-Body, Press-over-Metal

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GC UNITED KINGDOM Ltd.
22-23, Coopers Court
Newport Pagnell
UK - Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
E-mail info@uk.gceurope.com
www.uk.gceurope.com

GC FRANCE s.a.s.
9 bis, Avenue du Bouton d'Or - BP 166
F - 94384 Bonneuil sur Marne Cedex
Tel. +33.1.49.80.37.91
Fax. +33.1.49.80.37.90
E-mail info@france.gceurope.com
www.france.gceurope.com

GC EUROPE N.V.
Iberic Branch
Edificio Codesa 2
Playa de las Americas, 2, 1º, Of. 4
ES - 28230 Las Rozas, Madrid
Tel. +34.916.364.340
Fax. +34.916.364.341
E-mail info@spain.gceurope.com
www.spain.gceurope.com

GC EUROPE N.V.
Austrian Office
Tallak 124
A - 8103 Rein bei Graz
Tel. +43.3124.54020
Fax. +43.3124.54020.40
E-mail info@austria.gceurope.com
www.austria.gceurope.com

GC EUROPE N.V.
Benelux Office
Tooroplaan 11
NL - 3431 RC Nieuwegein
Tel. +31.30.604.88.87
Fax. +31.30.604.88.87
E-mail info@benelux.gceurope.com
www.benelux.gceurope.com

GC EUROPE N.V.
Finnish and Baltic States Office
Vanha Hommaksentie 11B
FIN - 02430 Masala
Tel. & Fax. +358.9.221.82.59
E-mail info@finland.gceurope.com
www.finland.gceurope.com

GC EUROPE N.V.
East European Office
Căzmanska 8
HR - 10000 Zagreb
Tel. +385.1.46.77.251
+385.1.46.78.474
Fax. +385.1.46.78.473
E-mail info@eeo.gceurope.com
www.eeo.gceurope.com

GC EUROPE N.V.
Scandinavian Office
Kungsporten 4 A
S - 427 50 Billdal
Tel. +46 31 939553
Fax. +46 31 914246
E-mail info@scandinavia.gceurope.com
www.scandinavia.gceurope.com

GC EUROPE N.V.
Swiss Office
Wilerstrasse 3
CH - 9545 Wängi
Tel. +41.52.366.46.46
Fax. +41.52.366.46.26
E-mail info@switzerland.gceurope.com
www.switzerland.gceurope.com

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Vita® is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.