

НОВИНКИ КАТАЛОГА 2020. БОРЫ МАНІ.

		Максимальный диаметр	Длина рабочая Длина общая	Макс. обороты
	Боры (5 шт) TC-72	ISO 165/010	7.6/21.6	300K
	Боры (5 шт) SF-S41	ISO109/010	3.6/16.6	300K
	Боры (5 шт) BC-42C	ISO 002/014	2.4/19.1	450K
	Боры (5 шт) BC-S42C	ISO 002/015	2.3/16.6	450K
	Боры (5 шт) BR-49C	ISO 001/08	0.8/19	450K
	Боры (5 шт) EX-11C	ISO 039/033	5/19	160K
	Боры (5 шт) FO-22C	ISO 298/017	7.1/21.1	450K
	Боры (5 шт) FO-S21C	ISO 257/015	7.1/18.2	300K
	Боры (5 шт) TC-21C	ISO 165/015	8.1/22	300K
	Боры (5 шт) TF-21C	ISO 171/017	7.1/21	450K
	Боры (5 шт) TF-S21C	ISO 171/017	7.1/18.8	450K
	Боры (5 шт) TO-58C	ISO 297/014	6/21.2	300K
	Боры (5 шт) TR-14C	ISO 198/023	9.2/22	300K
	Боры (5 шт) TR-21C	ISO 197/017	7.1/21.1	450K
	Боры (5 шт) TR-24C	ISO 197/019	7.1/21.2	300K
	Боры (5 шт) TR-256C	ISO 198/015	8.3/22.2	300K
	Боры (5 шт) EX-31F	ISO 234/013	4.1/19	300K
	Боры (5 шт) TR-12EF	ISO 199/014	10/21.7	300K
	Боры (5шт) CE-13	ISO 544/021	5.2/18.9	300K
	Боры (5шт) CE-17	ISO 172/012	8/21.9	300K