




DiaDent[®]

Serving 109 countries around the world ●●●

DiaDent Europe B.V.

©Copyright DiaDent Europe B.V.
www.diadenteurope.com



Welcome to DiaDent.

DiaDent is one of the leading Endodontic supply companies in the dental industry since 1985. DiaDent has been exporting different types of Gutta Percha Points and Paper Points to over 900 dealers in 109 countries all around the globe.

All of DiaDent's products are inspected by a strict Quality Management System following the ISO 9001 & ISO 13485 standards. With our promise and dedication to provide you with the finest dental products, you can always put your trust in DiaDent.

Always striving to be better, DiaDent puts in our utmost effort to improve our Endodontic lines such as the patented millimeter-marked Gutta Percha Points and Paper Points. Furthermore, we have expanded our product lines adding various new items to fulfill and satisfy our customers' needs and demands. DiaDent will continue to give our customers quality, value and services that far exceed their expectations.

Everyone here at DiaDent would like to thank you, our valued customers and business partners, for your continued support and loyalty. We wish everyone another happy and prosperous year ahead.

Sincerely,

DiaDent Group International



DiaDent Network

DiaDent Group International(Korea & International)

19-9 Songjung-Dong Heungduk-Gu Cheongju-si Choongchong Buk Do, Rep. of Korea
Tel)82-43-266-2315 Fax)82-43-262-8658
E-mail:diadent@chollian.net

DiaDent Group International Inc.(Canada & USA)

11-3871 North Fraser Way, Burnaby B.C, Canada V5J 5G6
Tel)1-604-451-8851 Fax)1-604-451-8865
E-mail:diadent@diadent.com

DiaDent Europe B.V (Europe)

Antennestraat 70, 1322 AS Almere, The Netherlands
Tel)31-36-549-8607 Fax)31-36-536-7317
E-mail:diadent@diadent europe.com

INDEX

ROOT CANAL FILLING MATERIALS • 6_11

- Millimeter Marked Gutta Percha Points
- Special Taper Gutta Percha Points
- Temporary Stopping
- Gutta Percha Gauge
- Millimeter Marked Paper Points
- Special Taper Paper Points
- Endo Organizers
- DiaPaste (Premixed Calcium Hydroxide Paste with Barium Sulfate)

ROOT CANAL OBTURATION • 12

- Wireless Obturation Device
- Gutta Percha Obturator

ENDODONTIC FILES • 13_16

- EDTA Root Canal Preparation Water
- EDTA Root Canal Preparation Cream
- Nickel Titanium Hand files
- Endo Stops & Dispenser
- Nickel Titanium Rotary Files
- Helix Organizer

SMALL EQUIPMENTS • 17_18

- Electronic Apex Locator
- LED Curing Light
- Intral Oral Camera

COSMETIC & RESTORATIVE • 19_22

- Restorative Hybrid Composite Resin
- Flowable Restorative Composite Resin
- Self Curing Composite Core Build-up
- Light-Cure Pit & Fissure Sealant
- Bonding Agent
- Etching Gel
- Light-Cure Temporary Filling Materials

INFECTION CONTROL • 23_24

- Protective Shields & Glass
- Sterilizing Containers and Trays

BUR BLOCKS & BUR STERILIZATION • 25_26

- Bur Blocks
- Magnetic Bur Blocks
- Bur Sterilization Boxes

PROPHY BRUSHES & CUPS • 27_28

- Prophy Brushes
- Prophy Cups

ENDODONTICS • 29_33

- Aluminum Endo Storage Boxes
- Endo Gauges
- Suction & Irrigation Syringe

INSTRUMENTS • 34_38

- Procedure Set up Trays & Tubes
- Instrument Cassettes
- Instrument Storage Boxes
- Wax Knives
- Triple Bite Tray

MISCELLANEOUS ITEMS • 39_43

- Dental Waxes
- Nickel Titanium Arch Wire
- Glass Dappen Dishes
- Medicament Bottle & Stand
- Glass Mixing Slab
- Retractor & Endo Training Bloc
- Amalgam Well

ROOT CANAL FILLING MATERIALS


DiaDent's Millimeter Marked Gutta Percha Points & Paper Points

- All DiaDent Points are Hand Jig Rolled to exacting standards.
- They are uniformly tapered and unsurpassed for consistency in sizing.
- All Points are laser inspected to ensure the tightest tolerances surpassing ISO standards.

• Millimeter Marked Gutta Percha Points


- Quick and easy method for depth measurements.
- Points with mm marking at 16mm, 18mm, 19mm, 20mm, 22mm and 24mm.
- Eliminates apex perforations which may lead to bleeding.
- Reaches the apex, completely seal the entire length of the canal.
- Save chair time with premeasured points.
- Most accurate depth measurements.

Millimeter Marked



• Computerized Laser Inspection System


- Each point is inspected by the laser beam.
- Using a laser beam, each DiaDent ML.029 point is individually inspected for diameters at D3 and D16.
- ISO and ADA standard requires the use of 1/100mm incremental measurement units, DiaDent uses 1/1,000mm increments.




• Millimeter Marked Paper Points

- Quick and easy method for depth measurements.
- Points with mm marking at 18mm, 19mm, 20mm, 22mm and 24mm.
- Eliminates apex perforations which may lead to bleeding.
- Completely drying the entire length of the canal.
- Save chair time with premeasured points.
- Most accurate depth measurements.

Millimeter Marked



DiaDent Quality Assurance




DEVIATION FROM STANDARD (mm)

■ DIADENT ML.029 SPECS
■ ISO and ADA SPECS (#10 – #25 sizes)
■ ISO and ADA SPECS (#30 – #140 sizes)

US Patent No. 6,004,133
International Patents pending.
European Patents No. 0964656
United Kingdom Patents No. 0964656
Deutsches Patents No. EP 0964656

Comparative absorption of Paper Points, by Brand




JOSE PUMAROLA-SUNE,
JOURNAL OF ENDODONTICS, VOL 24, NO.12, 796-798

Hand Jig Rolled Points vs Machine Rolled Points


DiaDent's method of hand rolling the points ensures the following:

DiaDent Hand Jig Rolled Gutta Percha Points



- Consistent, Uniform taper
- Smooth, Rounded Bullet shaped tips
- Readily conforms to shape of canal
- Greater likelihood of optimal apical seal

Machine Rolled Gutta Percha Points



- Inconsistent taper
- Uneven, Unpredictably shaped tips
- Overlapping of material of the tips
- Reduced likelihood of optimal apical seal

ROOT CANAL FILLING MATERIALS

GUTTA PERCHA POINTS

ML.029

MILLIMETER MARKED + LASER INSPECTED GUTTA PERCHA POINTS

ML.029

M - Millimeter Marked
L - Laser Inspected
.029 - Diameter Tolerance



• Package Information

• P-43 Vial Package

Contains 6 light weight plastic square vials with color coded caps
Box dimension : 77 x 17 x 50mm(W x H x D)

• P-63 Spill Proof Package

6 compartments with outer paper box
Box dimension: 73 x 17 x 49mm(W x H x D)

• Available Sizes

SIZE	P-43	P-63
#8	ML 101-401	ML 101-601
#10	ML 101-402	ML 101-602
#15	ML 101-403	ML 101-603
#20	ML 101-404	ML 101-604
#25	ML 101-405	ML 101-605
#30	ML 101-406	ML 101-606
#35	ML 101-407	ML 101-607
#40	ML 101-408	ML 101-608
#45	ML 101-409	ML 101-609
#50	ML 101-410	ML 101-610
#55	ML 101-411	ML 101-611
#60	ML 101-412	ML 101-612
#70	ML 101-413	ML 101-613
#80	ML 101-414	ML 101-614
#90	ML 101-415	ML 101-615
#100	ML 101-416	ML 101-616
#110	ML 101-417	ML 101-617
#120	ML 101-418	ML 101-618
#130	ML 101-419	ML 101-619
#140	ML 101-420	ML 101-620
#15/40	ML 101-491	ML 101-691
#45/80	ML 101-492	ML 101-692
#90/140	ML 101-493	ML 101-693

(Pkg. of 120)

SIZE	P-63
XF	ML 102-601
FF	ML 102-602
MF	ML 102-603
F	ML 102-604
FM	ML 102-605
M	ML 102-606
ML	ML 102-607
L	ML 102-608
XL	ML 102-609
XF/M	ML 102-691

(Pkg. of 100)

• Size vs Color Chart

15	45	90	XF	ML
20	50	100	FF	L
25	55	110	MF	XL
30	60	120	F	08
35	70	130	FM	10
40	80	140	M	

ROOT CANAL FILLING MATERIALS

SPECIAL TAPER GUTTA PERCHA POINTS

DIA-PROISO .04/.06 PLUS MILLIMETER MARKED + LASER INSPECTED SPECIAL TAPER GUTTA PERCHA POINTS

- Millimeter Marked Special Taper Gutta Percha Points.
- Designed to be used to complement for Profile ISO Rotary Files & K3 Files.
- Available in taper .04 & .06



Available Sizes

Size	Dia-ProISO .04 Plus	Dia-ProISO .06 Plus
	P-63	P-63
#15	ML 158-603	ML 159-603
#20	ML 158-604	ML 159-604
#25	ML 158-605	ML 159-605
#30	ML 158-606	ML 159-606
#35	ML 158-607	ML 159-607
#40	ML 158-608	ML 159-608
#45	ML 158-609	ML 159-609
#50	ML 158-610	ML 159-610
#55	ML 158-611	ML 159-611
#60	ML 158-612	ML 159-612
#70	ML 158-613	ML 159-613
#80	ML 158-614	ML 159-614
#15/40	ML 158-691	ML 159-691
#45/80	ML 158-692	ML 159-692

(Pkg. of 60)

DIA-PROT

- Designed to be used to complement **Protaper Files**.
- Millimeter Marked Points
- One hand slide open lid



Size	Item No.
DIA-PROT #F1	ML 150-S601
DIA-PROT #F2	ML 150-S602
DIA-PROT #F3	ML 150-S603
DIA-PROT #F4	ML 150-S604
DIA-PROT #F5	ML 150-S605
DIA-PROT #F1/F3	ML 150-S691
DIA-PROT #F4/F5	ML 150-S692
DIA-PROT #F1/F5	ML 150-S693

(Pkg. of 60)

DIA-PINK

- Shorter than usual to suit slightly shorter canals
- Available in 20mm length
- Widely used in Europe
- One hand slide open lid



Size	Item No.
#XXF	103-901
#XF	103-902
#F	103-903
#M	103-904
#L	103-905
#XXF/M	103-991

(Pkg. of 100)

FEATHERED TIP

- Designed to be used as auxiliary points to the master cone technique
- Extra fine taper
- Hand Jig Rolled



Size	Item No.
#XF	126-601
#FF	126-602
#MF	126-603
#F	126-604
#FM	126-605
#M	126-606
#ML	126-607
#L	126-608
#XL	126-609
#XF/M	126-691

(Pkg. of 100)

DIA-GTP

- Designed to be used to complement **Profile GT Greater Taper Files**.
- Available in .06, .08, .10, .12, taper.



Size	Item No.
DIA-GTP .06 #20	112-601
DIA-GTP .08 #20	112-602
DIA-GTP .10 #20	112-603
DIA-GTP .12 #20	112-604
DIA-GTP .06/.12 #20	112-691

(Pkg. of 60)

TEMPORARY STOPPING

- For temporary filling and thermal diagnostics tool.
- A warmed stick with pinpoint thermal sensitivity.
- High quality Gutta Percha.
- Easy to heat.
- Available in White, Yellow, Pink and Assorted.



	White	Yellow	Pink	Assorted
30g	180-101	180-102	180-103	180-191
120g	180-201	180-202	180-203	180-291
250g	180-301	180-302	180-303	180-391
500g	180-401	180-402	180-403	180-491



DIA-GP GAUGE

- ITEM #781-001
- Designed to size a perfect master cone from any Gutta Percha Points.
- Drop it into the apical size and cut off with a blade.
- Autoclavable
- Dimensions: 92 x 3 x 20mm/m(W x H x D)



ROOT CANAL FILLING MATERIALS

PAPER POINTS

MMPP

MILLIMETER MARKED PAPER POINTS

- MMPP is marked to provide preliminary guidance in measuring the root canal depth before proceeding clinically with the application of the Gutta Percha Points.
- These points are highly absorbent and color coded to identify ISO tip sizes.



· Package Information

· P-33 Cell Package

Points remain sterile until opened
Box dimension: 78 x 19 x 50mm (W x H x D)

· P-43 Vial Package

Contains 6 light weight plastic square vials
with color coded caps
Box dimension: 77 x 17 x 50mm (W x H x D)

· P-63 Spill Proof Package

6 compartments plastic package
Box dimension: 73 x 17 x 49mm (W x H x D)

· Available Sizes

SIZE	P-33	P-43	P-63
#8	MP 201-301	MP 201-401	MP 201-601
#10	MP 201-302	MP 201-402	MP 201-602
#15	MP 201-303	MP 201-403	MP 201-603
#20	MP 201-304	MP 201-404	MP 201-604
#25	MP 201-305	MP 201-405	MP 201-605
#30	MP 201-306	MP 201-406	MP 201-606
#35	MP 201-307	MP 201-407	MP 201-607
#40	MP 201-308	MP 201-408	MP 201-608
#45	MP 201-309	MP 201-409	MP 201-609
#50	MP 201-310	MP 201-410	MP 201-610
#55	MP 201-311	MP 201-411	MP 201-611
#60	MP 201-312	MP 201-412	MP 201-612
#70	MP 201-313	MP 201-413	MP 201-613
#80	MP 201-314	MP 201-414	MP 201-614
#90	MP 201-315	MP 201-415	MP 201-615
#100	MP 201-316	MP 201-416	MP 201-616
#110	MP 201-317	MP 201-417	MP 201-617
#120	MP 201-318	MP 201-418	MP 201-618
#130	MP 201-319	MP 201-419	MP 201-619
#140	MP 201-320	MP 201-420	MP 201-620
#15/40	MP 201-391	MP 201-491	MP 201-691
#45/80	MP 201-392	MP 201-492	MP 201-692
#90/140	MP 201-393	MP 201-493	MP 201-693

SIZE	P-33	P-63
XXF	MP 202-301	MP 202-601
XF	MP 202-302	MP 202-602
F	MP 202-303	MP 202-603
M	MP 202-304	MP 202-604
C	MP 202-305	MP 202-605
XC	MP 202-306	MP 202-606
XXF/XC	MP 202-391	MP 202-691

· Size vs Color Chart

8	10		
15	45	90	XXF
20	50	100	XF
25	55	110	F
30	60	120	M
35	70	130	C
40	80	140	XC

· Contents :

Size	P-33	P-43	P-63
#08-80	200pcs	200pcs	200pcs
#90-140	200pcs	120pcs	120pcs

ROOT CANAL FILLING MATERIALS

SPECIAL TAPER PAPER POINTS

DIA-PROISO .04/.06 PLUS

MILLIMETER MARKED SPECIAL TAPER PAPER POINTS



- Millimeter marked Special Taper Absorbency Paper Points.
- Designed to be used for drying and cleaning the canals prepared with Profile ISO Rotary Files & K3 Files.
- Available in .04 & .06 taper.

Size	Dia-ProISO .04 Plus	Dia-ProISO .06 Plus
	P-63	P-63
#15	MP 258-603	MP 259-603
#20	MP 258-604	MP 259-604
#25	MP 258-605	MP 259-605
#30	MP 258-606	MP 259-606
#35	MP 258-607	MP 259-607
#40	MP 258-608	MP 259-608
#45	MP 258-609	MP 259-609
#50	MP 258-610	MP 259-610
#55	MP 258-611	MP 259-611
#60	MP 258-612	MP 259-612
#70	MP 258-613	MP 259-613
#80	MP 258-614	MP 259-614
#15/40	MP 258-691	MP 259-691
#45/80	MP 258-692	MP 259-692

(Pkg. of 100)



DIA-PROT

- Designed to be used to complement Protaper Files.
- Millimeter Marked Points
- One hand slide open lid



Size	Item No.
DIA-PROT #F1	MP 250-S601
DIA-PROT #F2	MP 250-S602
DIA-PROT #F3	MP 250-S603
DIA-PROT #F4	MP 250-S604
DIA-PROT #F5	MP 250-S605
DIA-PROT #F1/F3	MP 250-S691
DIA-PROT #F4/F5	MP 250-S692
DIA-PROT #F1/F5	MP 250-S693

(Pkg. of 100)

FEATHERED TIP

- Designed to be used as auxiliary points to the master cone technique.
- Extra fine taper.
- Hand Jig Rolled.



Size	Item No.
#XXF	204-601
#XF	204-602
#F	204-603
#M	204-604
#C	204-605
#XC	204-606
#XXF/XC	204-691

(Pkg. of 200)

DIA-GTP

- Designed for drying and cleaning canals prepared with Profile GT Greater Taper Files.
- Available in .06, .08, .10, .12 taper.



Size	Item No.
DIA-GTP .06 #20(100/box)	208-601
DIA-GTP .08 #20(100/box)	208-602
DIA-GTP .10 #20(60/box)	208-603
DIA-GTP .12 #20(60/box)	208-604
DIA-GTP .06/.12 #20(60/box)	208-691

Package Information



P-33 Cell Package

Points remain sterile until opened
Box dimension:
78 x 19 x 50mm(W x H x D)



P-43 Vial Package

Contains 6 light weight plastic square vials with color coded caps
Box dimension:
77 x 17 x 50mm(W x H x D)



P-63 Spill Proof Package

6 compartments with outer paper box
Box dimension:
73 x 17 x 49mm(W x H x D)



SP-60 Slide Package

One hand slide open lid
Box dimension:
70 x 13 x 45 mm (W x H x D)

ROOT CANAL FILLING MATERIALS

ENDO ORGANIZERS

ENDO ORGANIZERS

- Packed in plastic Endo Organizer package with 12 compartments.
- The ideal way to organize Endodontic Points.
- DiaDent Endo Organizers can be ordered empty or a combination of Gutta Percha and Paper Points

ML.029 & MMPP

- Millimeter marked Gutta Percha and Paper Points.
- All DiaDent points are Hand Jig Rolled.



mmGP1



mmPP2



mmCOMBI1

Size	Description	Item No.
mmEG1	Empty box	EMG 301-100
mmGP1	240points per box(size15/40)	EMG 301-101
mmGP2	240points per box(size45/80)	EMG 301-102
mmGP3	240points per box(size15/80)	EMG 301-103
mmPP1	400points per box(size15/40)	EMP 301-201
mmPP2	400points per box(size45/80)	EMP 301-202
mmPP3	400points per box(size15/80)	EMP 301-203
mmCOMBI1	GP 120points/PP 200points per box(size15/40)	EMC 301-301
mmCOMBI2	GP 120points/PP 200points per box(size45/80)	EMC 301-302
mmCOMBI3	GP 240points/PP 400points per box(size15/40)	EMC 301-303
mmCOMBI4	GP 240points/PP 400points per box(size45/80)	EMC 301-304

.04 & .06 Tapered MMGP & MMPP


- Millimeter marked special taper Gutta Percha and Paper Points.
- Designed to be used to complement for Profile ISO Rotary Files & K3 Files.



mmGP2(6%)



mmPP1(4%)



mmCOMBI1(4%)

Size	Description	Item No.
mmGP 1(4%)	240points GP DIA-PROISO-.04 PLUS(size 15/80)	MSG 301-101
mmGP 2(6%)	240points GP DIA-PROISO-.06 PLUS(size 15/80)	MSG 301-102
mmGP 3(4% & 6%)	240points GP DIA-PROISO-.04 & .06 PLUS (size 15/40)	MSG 301-103
mmPP 1(4%)	240points PP DIA-PROISO-.04 PLUS(size 15/80)	MSP 301-201
mmPP 2(6%)	240points PP DIA-PROISO-.06 PLUS(size 15/80)	MSP 301-202
mmPP 3(4% & 6%)	240points PP DIA-PROISO-.04 & .06 PLUS (size 15/40)	MSP 301-203
mmCOMBI 1 (4%)	240points GP & PP DIA-PROISO-.04 PLUS(size15/40)	MSC 301-301
mmCOMBI 2 (6%)	240points GP & PP DIA-PROISO-.06 PLUS(size15/40)	MSC 301-302

ROOT CANAL FILLING MATERIALS

DIAPASTE™

CALCIUM HYDROXIDE PASTE WITH BARIUM SULFATE



Characteristics


- Premixed Calcium Hydroxide + Barium Sulfate paste for root canal treatment.
- Temporary root filling.
- Root canal filling for primary teeth.
- Use after pulp extirpation treatment.
- Filling on the pulpotomy.
- Temporary Pulp Capping.

Benefits

- Water Soluble: Easy to clean and remove, reduce operative time
- Easy and quick access
- Radiopaque
- Antibacterial

Packages

- ITEM #1001- 401 Regular Kit Type A - 2g of paste and 20 disposable tips with cap
- ITEM #1001- 402 Regular Kit Type B - 2g of paste and 10 disposable tips with cap
- ITEM #1001-315 Blue Disposable Syringe Tips – 20 disposable tips in a **blister package**



ROOT CANAL OBTURATION

Wireless Obturation Device

NEW

DIA-GUN

- Dia-Gun is Wireless Obturation System Device, Dia-Gun offers fast and easy obturation delivery. With this effective obturation system and technique, optimum results are achieved. The high-performance Dia-Gun device has the following features:

- Cordless and Compact Design allows complete freedom of movement
- Reliable Gutta Percha Filling to the desired level without delay
- Fast and continuous flow of filling tightly seals the canal
- Variable Temperature Settings (160°C, 180°C, 200°C): precise control of obturation flow
- Quick Heating System: temperature reaches 200°C within 25 seconds
- Ergonomically designed 360°C Swivel Needle offers improved access
- Thinner needle can easily fill the narrow root canal

Packages

- ITEM #4001-1001 Regular Kit
- Dia-Gun Handpiece & Charger
- Accessories: Safety Cap(2ea), Plunger, Plunger O-Rings(4ea), Power Adapter & Cord, Disposable Gun Tips(23G,4ea), Multipurpose wrench, Scrubbing Brush, Battery, Gutta Percha Obturator(100ea)

Accessories

- ITEM #4001-1010 Handpiece
- ITEM #4001-1020 Charger
- ITEM #4001-1030 Safety Cap(2ea)
- ITEM #4001-1040 Plunger
- ITEM #4001-1050 Plunger O-Rings(4ea)
- ITEM #4001-1060 Power Adapter & Cord
- ITEM #4001-1070 Disposable Gun Tips(23G,4ea)
- ITEM #4001-1071 Disposable Gun Tips(25G,4ea)
- ITEM #4001-1080 Multipurpose Wrench
- ITEM #4001-1090 Scrubbing Brush
- ITEM #4001-1100 Battery
- ITEM #181-101 Gutta Percha Obturator(Medium type, 100ea)



GUTTA PERCHA OBTURATOR

- ITEM #181-100 Soft type
- ITEM #181-101 Medium type
- ITEM #181-102 Hard type
- Designed to use with Dia-Gun the obturation gun.
- Specially formulated Gutta Percha with excellent flowing characteristics.
- Easy to heat
- Pkg. of 100



ENDODONTIC FILES

EDTA

Endo-Clear™

WATER BASED EDTA ROOT CANAL CLEANSER

Endo-Clear™ is the EDTA root canal cleanser, colorless and clear without foreign matters and effectively cleans the canal and removes smear layer. Chelating effect of Endo-clear™ softens root canal wall, dissolves calcium and removes it.

Characteristics

- Debridement in the root canal
- Dissolution of pulp remnants
- Removal of smear layer
- Microbicidal treatment within the canal
- Water based type

Benefits

- Chelating effect of Endo-Clear™ make root canal wall softened and calcium dissolved.
- Easy to enlarge the canal and clean it.
- Endo-Clear™ perfectly removes smear layer.
- Superior microbicidal treatment within the canal
- Excellent dissolution effect of pulp remnants
- Being based on water, more easier to get to apical area where hard for cream type to reach. So it effectively removes smear layer around apical area better than cream type.

Packages

Regular Package

- ITEM # 2002-5001 Regular Package EDTA 17% Solution 100 ml x 1
Dia-Flex (=Irrigation syringe) x 1
Disposable tips x 3
- ITEM # 2002-5002 Refill Package EDTA 17% Solution 100 ml x 1



DIA-PREP™ Plus

EDTA ROOT CANAL PREPARATION CREAM

EDTA Cream with Urea Peroxide, used as an irrigation agent for chemo-mechanical preparation of the root canal.

Characteristics

- Debridement in the root canal
- Dissolution of pulp remnants
- Removal of smear layer
- Facile movements of canal instruments
- Microbicidal treatment within the canal
- Clearance of dentin tubules and apical area

Benefits

- EDTA removes smear layer, softens dentin, and facilitates the removal of calcific obstruction.
- Urea Peroxide reacts with Sodium Hypochlorite to produce oxygen bubbles that flush out debris within the canal.
- The lubricant medium facilitates instrument movement within the canal to eliminate instrument binding or breakage.
- Unclogged dentin tubules and apical area permit better penetration of root canal medicaments and fillings.
- Provides rapid cleansing and enlarging of canals due to EDTA, which softens walls and smooth calcium deposits.
- Easier placement into canals than liquids.

Packages

- ITEM # 2002-1100 - 2 syringes of 6g, 10 disposable tips and 1 protective cap
- ITEM # 2001-6111 - 20 Disposable tips in a polybag



ENDODONTIC FILES

ENDODONTIC FILES

ENDODONTIC NICKEL TITANIUM HAND FILES / PLASTIC PACKAGE

- Increased Flexibility

: NiTi Shaft is made of premium Nickel-Titanium with 5 times the flexibility compared to Stainless Steel Files.

- Strength

: NiTi has 10 times the stress resistance compared to Stainless Steel Files.

- Convenience

: Produces smooth and progressive enlargement of root canal.

NiTi follow existing canal path and has memory so there is no need to Pre-Bend the files.

- Efficiency

: Color Coded handles and pre-inserted silicone stops.

- Safety

: NiTi resists spontaneous breakage.

Increased angular rotation reduces fracture.

NiTi decreases zipping and stripping.



NITI REAMERS



NICKEL-TITANIUM REAMERS(.02 TAPER)

SIZE	21MM	25MM
#8	NT501-101	NT501-201
#10	NT501-102	NT501-202
#15	NT501-103	NT501-203
#20	NT501-104	NT501-204
#25	NT501-105	NT501-205
#30	NT501-106	NT501-206
#35	NT501-107	NT501-207
#40	NT501-108	NT501-208
#45	NT501-109	NT501-209
#50	NT501-110	NT501-210
#55	NT501-111	NT501-211
#60	NT501-112	NT501-212
#70	NT501-113	NT501-213
#80	NT501-114	NT501-214
#15/40	NT501-191	NT501-291
#45/80	NT501-192	NT501-292



NITI K-FILES



NICKEL-TITANIUM K-FILES(.02 TAPER)

SIZE	21MM	25MM
#8	NT502-101	NT502-201
#10	NT502-102	NT502-202
#15	NT502-103	NT502-203
#20	NT502-104	NT502-204
#25	NT502-105	NT502-205
#30	NT502-106	NT502-206
#35	NT502-107	NT502-207
#40	NT502-108	NT502-208
#45	NT502-109	NT502-209
#50	NT502-110	NT502-210
#55	NT502-111	NT502-211
#60	NT502-112	NT502-212
#70	NT502-113	NT502-213
#80	NT502-114	NT502-214
#15/40	NT502-191	NT502-291
#45/80	NT502-192	NT502-292



NITI H-FILES



NICKEL-TITANIUM H-FILES(.02 TAPER)

SIZE	21MM	25MM
#8	NT503-101	NT503-201
#10	NT503-102	NT503-202
#15	NT503-103	NT503-203
#20	NT503-104	NT503-204
#25	NT503-105	NT503-205
#30	NT503-106	NT503-206
#35	NT503-107	NT503-207
#40	NT503-108	NT503-208
#45	NT503-109	NT503-209
#50	NT503-110	NT503-210
#55	NT503-111	NT503-211
#60	NT503-112	NT503-212
#70	NT503-113	NT503-213
#80	NT503-114	NT503-214
#15/40	NT503-191	NT503-291
#45/80	NT503-192	NT503-292

ENDODONTIC FILES

ENDODONTIC FILES

ENDODONTIC NICKEL TITANIUM HAND FILES / BLISTER PACKAGE

- Blister package for single use



NITI REAMERS



NICKEL-TITANIUM REAMERS(.02 TAPER)

SIZE	21MM	25MM
#8	NTB501-101	NTB501-201
#10	NTB501-102	NTB501-202
#15	NTB501-103	NTB501-203
#20	NTB501-104	NTB501-204
#25	NTB501-105	NTB501-205
#30	NTB501-106	NTB501-206
#35	NTB501-107	NTB501-207
#40	NTB501-108	NTB501-208
#45	NTB501-109	NTB501-209
#50	NTB501-110	NTB501-210
#55	NTB501-111	NTB501-211
#60	NTB501-112	NTB501-212
#70	NTB501-113	NTB501-213
#80	NTB501-114	NTB501-214
#15/40	NTB501-191	NTB501-291
#45/80	NTB501-192	NTB501-292

NITI K-FILES



NICKEL-TITANIUM K-FILES(.02 TAPER)

SIZE	21MM	25MM
#8	NTB502-101	NTB502-201
#10	NTB502-102	NTB502-202
#15	NTB502-103	NTB502-203
#20	NTB502-104	NTB502-204
#25	NTB502-105	NTB502-205
#30	NTB502-106	NTB502-206
#35	NTB502-107	NTB502-207
#40	NTB502-108	NTB502-208
#45	NTB502-109	NTB502-209
#50	NTB502-110	NTB502-210
#55	NTB502-111	NTB502-211
#60	NTB502-112	NTB502-212
#70	NTB502-113	NTB502-213
#80	NTB502-114	NTB502-214
#15/40	NTB502-191	NTB502-291
#45/80	NTB502-192	NTB502-292

NITI H-FILES



NICKEL-TITANIUM H-FILES(.02 TAPER)

SIZE	21MM	25MM
#8	NTB503-101	NTB503-201
#10	NTB503-102	NTB503-202
#15	NTB503-103	NTB503-203
#20	NTB503-104	NTB503-204
#25	NTB503-105	NTB503-205
#30	NTB503-106	NTB503-206
#35	NTB503-107	NTB503-207
#40	NTB503-108	NTB503-208
#45	NTB503-109	NTB503-209
#50	NTB503-110	NTB503-210
#55	NTB503-111	NTB503-211
#60	NTB503-112	NTB503-212
#70	NTB503-113	NTB503-213
#80	NTB503-114	NTB503-214
#15/40	NTB503-191	NTB503-291
#45/80	NTB503-192	NTB503-292

ENDO STOP DISPENSER

ENDO STOP DISPENSER

ITEM #701-000



- Contains 120 Assorted Round stops.
- Easy to use.
- Features ruler for measuring.
- Endo Stop Refillable.



ENDO STOPS

- Available Round Type.
- 100 Stops per poly bag.



White	701-101
Yellow	701-102
Red	701-103
Blue	701-104
Green	701-105
Black	701-106
Assorted	701-191

ENDODONTIC FILES

ENDODONTIC FILES

HELIX NICKEL TITANIUM ROTARY FILES

Benefits

- Made of Superior Nickel Titanium Alloy for Engine Driven use.
- Offers a smooth and progressive enlargement of the canal.
- Helix NiTi Rotary Files help to ease around curved canals and make a tapered preparation.
- Helix NiTi Rotary Files use ISO Color identification and come in 21mm and 25mm lengths.
- With amazing flexibility, easily move through canals, especially severely curved canals.
- High resistance to cyclic fatigue and file separation.
- Offers outstanding flexibility, cutting ability and reduce user fatigue.
- Ergonomic handle provides comfort and excellent tactile feedback.
- Recommended Speed: 150~350r/min



Size	21MM		25MM	
	.04	.06	.04	.06
#20	NT505-104	NT506-104	NT505-204	NT506-204
#25	NT505-105	NT506-105	NT505-205	NT506-205
#30	NT505-106	NT506-106	NT505-206	NT506-206
#35	NT505-107	NT506-107	NT505-207	NT506-207
#40	NT505-108	NT506-108	NT505-208	NT506-208
#20/35	NT505-193	NT506-193	NT505-293	NT506-293




.04 Taper

1 Line



.06 Taper

2 Lines



NEW

HELIX ORGANIZER

- ITEM #806-103 Red
- ITEM #806-104 Blue
- Ideal for Organizing Helix Files
- Holds 10 Helix Files
- Removable insert for Cleaning
- Dimensions: 50 x 50 x 45mm (W x H x D)




SMALL EQUIPMENTS

SMALL EQUIPMENTS

NAVIROOT PLUS - ELECTRONIC APEX LOCATOR

NaviRoot is a root canal length measuring device used in dental and endodontic treatment.

- **Accurate Measurement in Any Canal Condition**
 - Greater accuracy and predictability in wet canals
 - Automatic calibration for very accurate diagnosis in any canal condition
- **Larger Screen**
 - Large screen with digital and graphical feed back
 - Can be connected to personal computer and can be built in the unit chair
 - : The cable goes into RS232 port of a computer after setting up application software
- **Longer Battery Life**
 - Auto Power Off function(Timer 8 minutes)
 - Easy to read battery power indication
 - Use AA type alkaline battery(3ea included in a package)
 - Can be used rechargeable battery
- **Hi-Tech New Model for the Dentist**
 - Advanced function produces clinical accuracy rate and predictability
 - Hi-Technology design has great stability and ergonomic applicability
 - Self diagnosis for the system failure
- **Simple & Easy to use**
 - Holder & Touch Type Probes included
 - Foldable screen adjusts to desired viewing angle



Packages

- ITEM # 2120-2001 NaviRoot Plus Regular Kit Black
ITEM #2120-2002 NaviRoot Plus Regular Kit White
- Apex Locator
 - Accessories: File holder, Cable, Lip Clip, Y-shaped stick, Application software, RS232 Cable

Accessories

- ITEM #2120-1010 File Holder
- ITEM #2120-1020 Cable
- ITEM #2120-1030 Lip Clip
- ITEM #2120-1040 Y-shaped Stick
- ITEM #2120-1050 Application software
- ITEM #2120-1060 RS232 Cable



<NaviRoot Plus connected to the personal computer>

SMALL EQUIPMENTS

SMALL EQUIPMENTS

Q-LITE

LIGHTWEIGHT CORDLESS LED CURING LIGHT

DiaDent Q-lite is a lightweight cordless LED curing light unit that emits high intensity 1,100mW/cm² blue light in the wavelength range 430 - 490nm. It is designed for optimum polymerization of light-curing dental materials.

Benefits

- A lightweight and durable aluminum hand-piece that dissipates heat quickly.
- High intensity 1,100mW/cm² light output that offers fast polymerization.
- Wide wavelength range of 430 - 490nm that cover all composite.
- An efficient lithium-ion battery that does not build a memory effect.
- Long battery life with automatic power off.
- A detachable battery pack that allows for easy maintenance at a lower cost.
- Convenient preset curing programs, such as Two-Step and Soft-Start curing.
- 4 types of timer setting modes(5, 10, 15, 20 sec.)



Green(5sec) : Maximum light output for 5 seconds

Blue(10sec) : Maximum light output for 10 seconds

Orange(15sec) : An initial 70% output of maximum light intensity for 5 seconds, followed by 10 seconds output at full power(Two-Step curing program)

Red(20sec) : An initial 10% output of maximum light intensity that gradually increases to full power in 10 seconds, followed by 10 seconds output at full power(Soft-Start curing program)

Packages

- ITEM # 2101-1001 : Regular Kit (Light Guide 8 x 11mm)
- ITEM # 2101-1002 : Regular Kit (Light Guide 11 x 11mm)
- ITEM # 2101-1110 : 1 Hand-piece(120g)
- ITEM # 2101-1120 : 1 Charger Base
- ITEM # 2101-1121 : 1 Hand-piece support
- ITEM # 2101-1130 : 1 Light Guide Turbo(8 x 11mm)
- ITEM # 2101-1131 : 1 Light Guide Turbo(11 x 11mm)
- ITEM # 2101-1150 : 1 Lithium-ion Battery
- ITEM # 2101-1160 : 3 Light Protectors
- ITEM # 2101-1170 : 1 Power supply(Adapter)
- ITEM # 2101-1171 : 1 Power Cord

DiaCam USB 2.0 INTRAORAL CAMERA

Excellent resolution

DiaCam displays 1280 x 1024(1.3 megapixels) high-resolution digital images.

Easy installation

The USB cable goes into the camera hand-piece and into any USB 2.0 port of a computer after setting up the USB driver and application software(option).

Effective power control

Pushing the power on/off button on the hand-piece enables the camera to turn on/off power or LED(press 1.0 sec. for power on/off and 0.1 sec. or LED on/off). Power will automatically shut down after 5 minutes when the camera is not used.

Complete focus control

DiaCam has the focusing button on the hand-piece to adjust the focus of the lens.

Convenient capture function

DiaCam also has the capture button on the hand-piece to capture images. The captured images can be saved onto computer through the application software.

System requirements for computer


Pentium IV - 1 Ghz or higher /512 MB RAM / USB 2.0 port/64 MB Video memory / Windows 2000 or XP or Vista / DirectX 9.0 or higher

Packages

- ITEM # 2110-1002
USB type Intraoral camera set
- Camera Hand piece
- Accessories : 5m USB cable, camera holder, disposable sheaths(300pcs)
- Application software(bundle) & Device driver

Accessories

- ITEM #2110-1010 Sheath(Vinyl cap) 300pcs
- ITEM #2110-1020 USB cable 2.0m
- ITEM #2110-1025 USB cable 5.0m



RESTORATIVE MICRO HYBRID COMPOSITE RESIN

DIAFIL™

A LIGHT CURE, ESTHETIC RESTORATIVE MICRO HYBRID COMPOSITE RESIN

Micro hybrid composite resin for use in Anterior and Posterior restorations

Benefits

- Unique fillers offer minimum polymerization shrinkage which helps reduce stress and sensitivity on the tooth potentially
- Has excellent fracture toughness, high tensile, compressive strengths for long lasting results
- Very easy handling of acceptable viscosity and good marginal adaptation
- Remarkable color stability and superior quality of esthetics
- High Level of radiopacity

Packages

· Start Kit


- Shade (5 syringes x 4g)
A set (A1, A2, A3, A3.5, B2)
B set (A2, A3, A3.5, B2, B3)
C set (A1, A2, A3, A3.5, C2)
- One bottle bonding agent(5ml)
- Etching gel 1 syringe(3ml)
- Accessories : 50 Brush tips, 1 brush handle, mixing well, mixing pad, 10 disposable needle tips

START KIT A SET	2001-1101
START KIT B SET	2001-1102
START KIT C SET	2001-1103



Properties of Composite Resin

Physical properties



· Refill Package

- 1 syringe x 4g


A1	2001-1201
A2	2001-1202
A3	2001-1203
A3.5	2001-1204
A4	2001-1205
B1	2001-1206
B2	2001-1207
B3	2001-1208
C2	2001-1210
C3	2001-1211
A20	2001-1212
A30	2001-1213

· Capsule Type

- 0.25g capsule x 20ea

A1	2001-1405
A2	2001-1401
A3	2001-1402
A3.5	2001-1403
B2	2001-1404

SEM Micrographs of Composite Fillers



A(USA)



B(USA)



DiaFil

DIAFIL™ FLOW

LIGHT CURED FLOWABLE RESTORATIVE COMPOSITE RESIN

Flowable composite with very easy handling and high strength

Benefits

- Has a high level of fluoride release with no slump, run which helps a fast techniques and handling
- Enables to be used in various restorative procedures
 - : Class III, V and IV direct restorations
 - A liner for Class I, II
 - Pit & Fissure sealant where is required
- Well suited for minimally invasive dentistry
- Excellent tensile strength and low shrinkage
- High radiopaque materials

Packages

· Economic Package

- Flowable Resin 2g x 4 syringes
(Available shade A1, A2, A3, A3.5, B2, B3, C2, A2O, A3O)
- 40 disposable tips(black)

· Refill Package


- Flowable Resin 2g x 1 syringe
(Available shade A1, A2, A3, A3.5, B2, B3, C2, A2O, A3O)
- 10 disposable tips(black)

SHADE	Economic Pkg.	Refill Pkg.
A1	2001-4101	2001-4201
A2	2001-4102	2001-4202
A3	2001-4103	2001-4203
A3.5	2001-4104	2001-4204
B2	2001-4105	2001-4205
B3	2001-4106	2001-4206
C2	2001-4107	2001-4207
A2O	2001-4108	2001-4208
A3O	2001-4109	2001-4209




Properties


Compressive Strength




Diametral tensile Strength



Flexural Strength



Flexural Modulus



DIAFIL CORE™

SELF-CURING COMPOSITE CORE BUILD-UP

Self-curing and long-term fluoride releasing composite materials consist of two-components, base plus catalyst for core build-up system

Benefits

- Provides good handling properties for easier use
- Offers more faster self cure set time
- Smooth and creamy materials help fill it into needle tips very easily
- High level of radiopacity

Packages

- ITEM #2001-7101
 - 10g Base + 10g Catalyst
 - 30 Spatulas
 - 50 Mixing pads
 - 10 Needle tips



DIASEAL™

A LIGHT CURE PIT & FISSURE SEALANT

Designed for sealing enamel pits and fissures of teeth to prevent caries

Benefits

- UDMA series monomer and Nano-filler offers high flowability, excellent wear resistance and strength
- Optimal viscosity to flow easily and helps for easy application
- Ensures hygienic and easy delivery and outstanding versatility
- High visual verification during application

Packages

Regular Package

ITEM #2001-5101

- Sealant 2 syringes x 1.2ml
- Etching 2 syringes x 1.2ml
- 10 disposable needle tips

Refill Package

ITEM #2001-5201

- Sealant 1 syringe x 1.2ml
- 5 disposable needle tips



Properties

Mean value and standard deviation of micro-leakage score

Materials	Micro-leakage score (standard deviation)
Dia-Seal™	1.0000(1.1094)
U	1.7333(1.4829)
M	0.9167(1.2795)

Thermocycling 500 times,
1% Methylene blue solution

Shear bond strength(MPa)



DIAPLUS™

SINGLE COMPONENT BONDING AGENT

Recommended for all direct composite restorations, light-cured bonding procedures

Benefits

- High dentin, enamel bond strength performance
- Available in all classes of direct composite restoration
- Indirect restorations of Porcelain, Amalgam or Composite
- Repair after light curing with a minimal film thickness
- Ethanol-based bonding agent enables to maintain its exceptional bond strength

Packages

- ITEM #2001-2101
- One-bottle bonding agent(1 bottle x 5ml)



COSMETIC & RESTORATIVE

PIT & FISSURE SEALANT & BONDING AGENT

DIAETCH™ 37%

A 37% ETCHING GEL PHOSPHORIC ACID

A phosphoric acid is used in etch enamel and/or dentin

Benefits

- Phosphoric acid etching of enamel and dentin with ideal thickness
- Quick and effective application to the desired place
- Enables to rinse completely and easily from the tooth surface
- High biocompatibility

Packages

• Economic Package Type A

ITEM #2001-3101

- 5 Syringes x 5ml and 20 disposable tips

• Economic Package Type B

ITEM #2001-3102

- 5 syringes x 3ml and 20 disposable tips

• Regular Package

ITEM #2001-3103

- 2 syringes x 3ml and 10 disposable tips

• Refill Package Type A

ITEM #2001-3202

- 1 syringe x 5ml and 5 disposable tips



Before

After

(Microphotographs of premolar, x1,000)

DIA-TEMP™

A VISIBLE, LIGHT CURED TEMPORARY FILLING MATERIAL

Light Cured material with nano-sized silver for temporary/intermediate restorations of inlay and onlay, for relining prefabricated crowns and bridges.

Benefits

- Quick and simple removal: Can be removed in one piece
- High wear resistance, mechanical strength
- Ensures a fast set
- Excellent adhesion
- Easy application: Convenient syringe dispensed material does not required mixing
- Blue Dia-Temp is ideal for Posterior teeth
- Yellow Dia-Temp is ideal for Anterior teeth

Packages

3 syringes of 3g

ITEM #2003-1101 Blue

ITEM #2003-1201 Yellow



INFECTION CONTROL

PROTECTIVE SHIELDS & GLASSES

Dia-100FS/ Dia-100SS

- Optical grade protective shields and glasses for face and eyes from flying debris or substances.
- Lightweight and very comfortable.
- Shields can be quickly positioned and Fits over eyeglasses.
- Easy to clean, washable.

• Dia-100FS

- Full face protection shield
- 1 Frame w/5 shields

1005-5010	DIA-100FS	REFILL PACKAGE	
1005-5011	DIA-100FS	REGULAR PACKAGE	CRYSTAL
1005-5012	DIA-100FS	REGULAR PACKAGE	WHITE
1005-5014	DIA-100FS	REGULAR PACKAGE	BLUE

• Dia-100SS

- Eye protection glasses
- 1 Frame w/5 shields

1005-5020	DIA-100SS	REFILL PACKAGE	
1005-5021	DIA-100SS	REGULAR PACKAGE	CRYSTAL
1005-5022	DIA-100SS	REGULAR PACKAGE	WHITE
1005-5024	DIA-100SS	REGULAR PACKAGE	BLUE

Anti-fog



Dia-100FS



Dia-100SS

Dia-300i

- ITEM #1005-1041
- Protection Glasses. Protect your eyes from flying debris during clinical or laboratory procedures.
- High optical quality lenses with chemical resistance properties.
- Comfort fit at the nose bridge at the temple.
- Lens coating(PC lenses, Anti-fog, Anti-UV, Anti-static)



Anti-fog

Anti-UV


Anti-static



Product Features



Vertical adjustments



Horizontal adjustments



Ultra light weight

Dia-200B

- ITEM #1005-1031
- Multi Purpose Goggles. Used among dentists and surgeons for eyes protection from flying debris. Preventing eye injuries and contamination from fine dust.
- Spherical lens design for superior optics.
- Comfortable lightweight fit.
- Fits over Eyeglasses.
- Designed for coolness and ventilation.
- Anti-fog coating, Anti-static.



Anti-fog

Anti-static



Product Features



Superior cushioning and ventilation



Good air circulation



Easily used with correction glasses

Dia-400D

- ITEM #1005-1051
- The protective glasses for blue curing light(400-510nm) are used to protect your eyes when curing composite resin.
- Excessive exposure to blue lights will damage your vision and your eyes health.
- Lenses color: Amber




Anti-fog

Anti-UV

Anti-scratch

Anti-static



Product Features



Vertical adjustments



Horizontal adjustments



Ultra light weight

STERILIZING CONTAINERS & TRAYS

GERMICIDE TRAY

- ITEM #711-001



- Used for Cold Sterilization.
- Easy & Safe to handle.
- One Large & One Small compartment.
- Easy loading and unloading.
- Dimensions: 263 x 72 x 115 m/m (W x H x D)



GERMICIDE TRAY TYPE B

- ITEM #711-102 Grey
- ITEM #711-104 Mauve
- ITEM #711-106 Blue



- Efficient high-volume germicide tray designed for maximum use of chemicals to clean/disinfect/sterilize instruments
- Dimensions : 350 x 150 x 250 m/m (W x H x D)



CONTAINER (STERILIZING)



- ITEM #824-001 Type A
 - Holds 10 Endo Instruments.
 - Dimensions : 185 x 35 x 145 m/m (W x H x D)
- ITEM #824-002 Type B
 - Holds 16 Endo Instruments.
 - Dimensions : 185 x 35 x 286 m/m (W x H x D)
- Well Built, Rigid, Sterilizing Container.
- Features Locking System, Fits many sterilizing units.



Type A



Type B

CONTAINER(STORAGE)



- ITEM #829-001 Type A
 - Dimensions: 255 x 50 x 216m/m (W x H x D)
- ITEM #829-002 Type B
 - Dimensions: 255 x 50 x 108m/m (W x H x D)
- Multi Purpose Box.
- Use as an organizer or for storage.
- Bur Holders fit well in this container.



Type A



Type B

STERILIZING DRUM



- ITEM # 828-001 Small
 - Diameter: 75m/m, Height: 80m/m
- ITEM # 828-002 Large
 - Diameter: 115m/m, Height: 100m/m

Autoclaving burs & smaller Instruments Drainage



Large



Small

WASTE CONTAINER



- ITEM #814-001 Small
 - Diameter: 56m/m, Height: 47m/m
 - Color: Silver
- ITEM #815-002 Large
 - Diameter: 70m/m, Height: 74m/m
 - Color: Gold

Separates into 3 pieces for easy Cleaning. Equipped with locking system.



Small




Large

FLAT TRAY(STERILIZING) TYPE A

- ITEM #827-001 Blue



- Sterilizing Tray.
- Ideal for Operator Instruments.
- Anodized Aluminum.
- Dimensions: 200 x 20 x 105m/m (W x H x D)



FLAT TRAY(STERILIZING) TYPE B

- ITEM #827-002 Red



- Sterilizing Tray.
- Ideal for Operator Instruments.
- Anodized Aluminum.
- Dimensions: 180 x 10 x 100m/m (W x H x D)



BUR BLOCKS & BUR STERILIZATION

BUR BLOCKS

BUR BLOCK W/COVER



Type A

· ITEM # 836-001



- Holds 24 Friction Grip (high speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 44m/m (W x H x D)



Type B

· ITEM # 836-002



- Holds 24 Right Angle (low speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 44m/m (W x H x D)



Type C

· ITEM # 836-003



- Holds 12 Right Angle (low speed) Burs & 12 Friction Grip (high speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 44m/m (W x H x D)



Type D

· ITEM # 836-004



- Holds 36 Friction Grip (high speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 50m/m (W x H x D)



Type E

· ITEM # 836-005



- Holds 36 Right Angle (low speed) Burs.
- Durable Aluminum with Plastic Cover.
- Remove Plastic Cover before Autoclave.
- Dimensions: 60 x 45 x 50m/m (W x H x D)

STERI BUR BLOCK W/COVER



Type A

· ITEM #837-001



- Holds 24 Friction Grip (high speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 44m/m (W x H x D)



Type B

· ITEM #837-002



- Holds 24 Right Angle (low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 44m/m (W x H x D)



Type C

· ITEM #837-003



- Holds 12 Friction Grip (high speed) Burs & 12 Right Angle (low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 44m/m (W x H x D)



Type D

· ITEM #837-004



- Holds 36 Friction Grip (high speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 50m/m (W x H x D)



Type E

· ITEM #837-005



- Holds 36 Right Angle (low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 x 45 x 50m/m (W x H x D)

BUR BLOCK

Type A

· ITEM # 832-001



- Holds 48 Friction Grip (high speed) Burs.
- Durable Aluminum.
- Dimensions: 77 x 7 x 40m/m (W x H x D)



Type B

· ITEM # 832-002



- Holds 48 Right Angle (low speed) Burs.
- Durable Aluminum.
- Dimensions: 77 x 7 x 40m/m (W x H x D)



Type C

· ITEM # 832-003



- Holds 24 Friction Grip (high speed) Burs & 24 Right Angle (low speed) Burs.
- Durable Aluminum.
- Dimensions: 77 x 7 x 40m/m (W x H x D)




BUR BLOCKS & BUR STERILIZATION

MAGNETIC BUR BLOCKS & BUR STERILIZATION

STERI BUR GUARD(22 HOLE)

- ITEM #716-103 Yellow
- ITEM #716-106 Blue
- Fits 22 latch & friction grip burs
- Lid does not allow the burs to fall out of the stand.
- Diameter: 75m/m, Height: 35m/m




MAGNETIC BUR BLOCK(SQUARE-8 HOLE)

- ITEM #717-303 Yellow
- ITEM #717-306 Blue
- Fits both FG and RA types.
- Dimensions : 40 x 13 x 40m/m (W x H x D)



MAGNETIC BUR BLOCK(OBLONG-14 HOLE)

- Type A
- ITEM #717-103 Yellow
 - ITEM #717-105 Beige
 - ITEM #717-106 Blue
 - Dimensions: 77 x 13 x 33m/m (W x H x D)



Type A

- Type B
- ITEM #717-203 Yellow
 - ITEM #717-205 Beige
 - ITEM #717-206 Blue
 - Dimensions: 76 x 12 x 37mm (W x H x D)
 - 14 Holes
 - Fits both FG and RA types.



Type B

MAGNETIC BUR BLOCK(ROUND 24 HOLE)




- 8 Holes
- ITEM #717-503 Yellow
 - ITEM #717-505 Beige
 - ITEM #717-506 Blue
 - Diameter : 43m/m, Height: 10m/m



8 Holes

- 24 Holes
- ITEM #717-603 Yellow
 - ITEM #717-605 Beige
 - ITEM #717-606 Blue
 - Diameter : 73m/m, Height: 10m/m
 - Fits both FG & RA types



24 Holes

GERMICIDE ROUND BUR BOX

- ITEM #915-001
- Contains perforated insert for drainage.
- Practical & Durable.
- 1 Basin bowl/ 1 Perforated insert/ 1 Cover.
- Diameter: 75m/m, Height: 35m/m



BUR STERI BATH

- ITEM # 730-001 White
- ITEM # 730-002 Transparent




- Used for Sterilizing Burs & Small Instruments.
- Removable insert for liquid drainage.
- 1 Basin bowl/ 1 Perforated insert/ 1 Cover.
- Diameter : 88 m/m
- Height : 68 m/m



AL BOX FOR BUR HOLDER(W/7 INSERTS)

- ITEM # 831-001
- Holds 24 High Speed Burs & 32 Low Speed Burs.
- Interchangeable Inserts. Versatile.
- Easy to transport.
- Dimensions : 125 x 53 x 130m/m (W x H x D)



DOWEL PIN

- ITEM #907-001: Twin pin
- ITEM #907-002: Single pin



Twin pin with metal sleeve & Rubber cap

- Rubber cap saves working time.
- Rubber cap facilitates work.
- Exact fit into preparations and no distortion.
- Pkg.of 1,000



PROPHY BRUSHES & CUPS

PROPHY BRUSHES

SCALING BRUSH LATCH TYPE

- ITEM #941-001
- Packed 72 brushes.
- Regular cup style.
- For polishing of tooth surface after scaling.
- Useful for wide space.
- Diameter: 4.49mm



CANDLE BRUSH

- ITEM #941-104 : Latch Type, Packed 24 brushes
- ITEM #941-105 : Screw Type, Packed 24 brushes and 10 mandrels
- For removing stains from occlusal pits and fissures and for preparing tooth surface prior to sealant and orthodontic application.
- Diameter
- Latch Type : 4.50mm
- Screw Type : 5.00mm



Candle Brush (Latch Type)



Candle Brush (Screw Type)

PROPHY BRUSH

- ITEM #950-101 : Latch Type
- ITEM #950-102 : Screw Type
- Packed 144 brushes
- Diameter
- Latch Type : 4.48mm
- Screw Type : 4.48mm



Latch Type (PB - 001)



Screw Type (PB - 002)

SCALING POLISHING BRUSH LATCH TYPE

- ITEM #941-103
- Packed 36 brushes.
- Small cup style.
- Useful for narrow space.
- Diameter : 2.27mm



POST BRUSH

- ITEM #941-101 : Latch Type, Packed 30 brushes
- ITEM #941-102 : Screw Type, Packed 36 brushes and 10 mandrels
- For quick cleaning root canal preparation, under bridge.
- Available to insert sterilizing liquid during turn it left in Contra Angle.
- Available to flow out the stains like cement, liquid during turn it right.
- Recommended rotation: 300-800rpm.
- Single use.
- Diameter
- Latch Type : 1.20mm
- Screw Type : 1.03mm



Post Brush (Latch Type)



Post Brush (Screw Type)

PROPHY BRUSH POINTED HEAD

- ITEM #950-111 : Latch Type
- ITEM #950-112 : Screw Type
- Packed 144 brushes
- Diameter
- Latch Type : 4.48mm
- Screw Type : 4.48mm



Latch Type (PB - 003)



Screw Type (PB - 004)

PROPHY BRUSH LARGE HEAD LATCH TYPE

- ITEM #950-121
- Packed 144 brushes
- Diameter : 5.98mm



(PB - 005)

PROPHY BRUSH LARGE POINTED HEAD LATCH TYPE

- ITEM #950-131
- Packed 144 brushes
- Diameter : 5.98mm



(PB - 006)

PROPHY BRUSHES & CUPS

PROPHY CUPS

SILICONE CUP FOR POLISHING LATCH TYPE

- ITEM #941-002
- Soft, white cups
- For polishing the tooth surface, achieves a high shine quickly.
- Use conveniently after insert in Contra Angle
- Pkg. of 72



PROPHY RUBBER CUP RIBBED & WEBBED

- ITEM #950-201 : Latch Type
- ITEM #950-202 : Screw Type
- Packed 144 cups



Latch Type (RC - 001)



Screw Type (RC - 002)

PROPHY RUBBER CUP WEBBED

- ITEM #950-211 : Latch Type
- ITEM #950-212 : Screw Type
- Packed 144 cups



Latch Type (RC - 003)



Screw Type (RC - 004)

PROPHY RUBBER CUP RIBBED

- ITEM #950-221 : Latch Type
- ITEM #950-222 : Screw Type
- Packed 144 cups



Latch Type (RC - 005)



Screw Type (RC - 006)

PROPHY RUBBER CUP UMBRELLA SHAPE

- ITEM #950-231 : Latch Type
- ITEM #950-232 : Screw Type
- Packed 144 cups



Latch Type (RC - 007)



Screw Type (RC - 008)

PROPHY RUBBER CUP WEBBED WITH WASHER LATCH TYPE

- ITEM #950-241
- Packed 144 cups



(RC - 009)

PROPHY RUBBER CUP SNAP-ON

- ITEM #950-251 : Ribbed
- ITEM #950-252 : Webbed
- Packed 144 cups



Ribbed (RC - 010)



Webbed (RC - 011)

PROPHY RUBBER CUP SNAP-ON RIBBED & WEBBED

- ITEM #950-253
- Packed 144 cups



(RC - 012)

ENDODONTICS

ALUMINUM ENDO STORAGE BOXES

SMALL ROUND TYPE A

- ITEM # 803-003 Red
- ITEM # 803-004 Blue



- Has 6 Silicone Color Coded Silicone Cups.
- Ideal for Organizing Gutta Percha Points & Paper Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Diameter: 40m/m, Height: 60m/m



SMALL ROUND TYPE B

- ITEM # 802-003 Red
- ITEM # 802-004 Blue



- Ideal for Organizing Hand Files.
- Holds 20 instruments.
- Removable Insert for Cleaning.
- Diameter: 40m/m, Height: 60m/m



SMALL SQUARE TYPE A

- ITEM #807-003 Red
- ITEM #807-004 Blue



- Has 6 Silicone Color Coded Silicone Cups.
- Ideal for Organizing Gutta Percha Points & Paper Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 50 x 50 x 45m/m (W x H x D)



SMALL SQUARE TYPE B

- ITEM #806-003 Red
- ITEM #806-004 Blue



- Ideal for Organizing Hand Files.
- Holds 30 Endo Instruments.
- Removable Insert for Cleaning.
- Dimensions: 50 x 50 x 45m/m (W x H x D)



MEDIUM OBLONG TYPE A

- ITEM #801-003 Red
- ITEM #801-004 Blue



- A Versatile & Practical Organizer with a Clean Grip for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Holds 36 Instruments, 3 Silicone Cups, 1 Clean Grip.
- Removable Insert for Cleaning.
- Dimensions: 130 x 50 x 50m/m (W x H x D)



MEDIUM OBLONG TYPE B

- ITEM #809-003 Red
- ITEM #809-004 Blue



- A Versatile & Practical Organizer for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Holds 60 Instruments, 4 Silicone Cups.
- Removable Insert for Cleaning.
- Dimensions: 130 x 50 x 50m/m (W x H x D)



MEDIUM ROUND

- ITEM # 821-001 Silver



- Ideal for Organizing Hand Files & Gutta Percha Points & Paper Points.
- Holds 64 Instruments, 4 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Diameter: 92m/m, Height: 69m/m



MEDIUM SQUARE

- ITEM #804-003 Red
- ITEM #804-004 Blue



- A Versatile & Practical Organizer with a Clean Grip for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Features a Endo Meter for measuring Files & Points.
- Holds 5 Burs, 30 Instruments, 6 Silicone Cups, 1 Clean Grip, Endo Meter and 1 location for small instruments.
- Removable Insert for Cleaning.
- Dimensions: 120 x 50 x 80m/m (W x H x D)



ENDODONTICS

ALUMINUM ENDO STORAGE BOXES

DIA-DELUXE

- ITEM #811-003 Red



- DiaDent's Ultimate Box for Organizing your Files & Points. Holds 4 Spreaders/ Pluggers, 56(21mm) & 48(25mm) Instruments, 6 Gates Glidens, 10 Barbed Broaches, 28 Cups, Section for Cotton.
- Removable Insert for Cleaning.
- Dimensions: 215 x 50 x 100m/m (W x H x D)



SMALL OBLONG- 6 COMPT.

- ITEM #805-003 Red
- ITEM #805-004 Blue



- Ideal for Holding Burs & Small Instruments.
- Has 6 Compartments.
- Removable Insert for Cleaning.
- Dimensions: 130 x 15 x 50m/m (W x H x D)




LARGE OBLONG TYPE A

- ITEM #812-004 Blue



- A Simple & Effective Organizer with Sections for your Hand Files.
- Holds 180 Endo Instruments, 8 Color Coded Silicone Cups for Points, perforated bottom case.
- Removable Insert for Cleaning.
- Dimensions: 218 x 53 x 105m/m (W x H x D)



LARGE OBLONG TYPE B

- ITEM #812-005



- An Organizer for Hand & Rotary Files as well as Points.
- Store your Instruments in One Place.
- Holds 357 Endo Instruments, 21 Silicone Cups for Points.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE OBLONG TYPE C

- ITEM #812-006



- A Very Practical Box. Easy to use inserts individually as needed.
- Holds 206 Endo Instruments, 12 Silicone Cups for Points.
- Interchangeable Inserts.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE OBLONG TYPE D

- ITEM #812-007



- Organizer for Hand & Rotary Files.
- Holds 192 Hand Instruments, 80 Rotary Instruments.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE OBLONG TYPE E

- ITEM #812-008



- Store your files in one place. Holds 308 Endo Instruments, 5 Color Coded Silicone Cups, 3 Locations.
- 2 Separate Inserts. Height Adjustment(not visible).
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE COMBI

- ITEM # 818-001



- Versatile & Practical. Allows separation of specific instruments.
- Holds 5 individual Inserts(A+B+C+D+E).
- Insert A : 6 cups + 6 cups
- Insert B : 2 sets of 6 cups & Holds 12 Endo Instruments
- Insert C : 4 cups, Holds 48 Endo Instruments
- Insert D : 2 cups, Holds 48 Endo Instruments
- Insert E : Holds 72 Endo Instruments
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



ENDODONTICS

ALUMINUM ENDO STORAGE BOXES

MEDIUM COMBI

· ITEM #817-001



- A Compact & Practical File & Point Organizer.
- Holds 3 Individual Inserts (C+D+E).
- Insert C : 4 cups, Holds 48 Endo Instruments
- Insert D : 2 cups, Holds 48 Endo Instruments
- Insert E : Holds 72 Endo Instruments
- Removable Insert for Cleaning.
- Dimensions: 122 x 50 x 100m/m (W x H x D)



SMALL COMBI TYPE A

· ITEM # 816-010



- A Compact organizer for Gutta Percha & Paper Points.
- Holds 12 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE B

· ITEM #816-020



- A Compact Organizer for Hand & Rotary Files.
- Hold 24 Endo Instruments, 12 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE C

· ITEM #816-030



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments, 4 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE D

· ITEM #816-040



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments, 2 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE E

· ITEM #816-050



- A Compact Organizer for Hand & Rotary Files.
- Holds 72 Endo Instruments.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE F

· ITEM #816-060



- A Compact Organizer for Gutta Percha & Paper Points.
- Holds 12 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE G

· ITEM #816-070



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Medium Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



ENDODONTICS

ALUMINUM ENDO STORAGE BOXES

SMALL COMBI TYPE H

- ITEM #816-080



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Large Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE I

- ITEM #816-090



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Medium & Large Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



ENDO GRIP(MINI)

- ITEM #833-001



- Commonly Used for Cleaning Files.
- Holds sponge or gauze.
- Available in Blue.
- Diameter: 40m/m, Height: 35m/m



ENDO CLEANER

- ITEM #834-001 Open
- ITEM #834-002 Closed



- Commonly Used for Cleaning and sterilizing Files.
- Top unscrews for easy replacement of Gauze.
- Diameter: 40m/m, Height: 30m/m




LARGE OBLONG EMPTY BOX

- ITEM #820-003 Blue



- An all purpose storage container.
- Ideal for holding accessories.
- Dimensions: 210 x 25 x 50m/m (W x H x D)




CANAL INSTRUMENT HOLDER

- ITEM #830-001



- Holds 18 Endo Instruments (Files, Reamers).
- 45 degree perforations.
- Dimensions: 50 x 25 x 50m/m (W x H x D)



ENDO METER SQUARE

- ITEM # 835-001



- Holds & Measures 32 Files or Other Endo Instruments.
- Gauge from 15 to 30.5m/m, including half measurement.
- 2 Cups for Files, 1 location for Endo Stops.
- Dimensions: 125 x 40 x 40m/m (W x H x D)



ENDO METER

- ITEM #835-002 Type A- Silver
- ITEM #835-003 Type B- Red



- Holds & Measures 10 Files.
- Gauges from 1 to 35mm.
- Features a locking spring mechanism inside for easy use.
- Diameter: 35m/m, Height: 48m/m



Type A

Type B

ENDO GAUGES & IRRIGATION SYRINGES

FILE MEASURE SET(PLASTIC)

- ITEM #780-001
 - Holds 32 Measuring Grades.
 - Gauge from 12 to 27.5m/m, Including half measurements.
 - Dimensions: 87 x 31 x 20m/m (W x H x D)



STERI ENDO GUARD W/RULER(ENGINE FILE)

- ITEM #716-203 Yellow
- ITEM #716-206 Blue
- Holds 16 NiTi Rotary Files.
- Right or left hand use.
- Available to insert sponge to clean or hold files.
- Diameter: 75m/m, Height: 35m/m



ENDO FINGERING

- ITEM #790-101
 - Provides complete control of entire endodontic files.
 - Disposable sponges hold files in a convenient position.
 - Autoclavable except cups
- Contents :
 - ITEM #790-111 : Disposable Cups(48pcs),
 - ITEM #790-121 : Disposable Sponges(48pcs),
 - Finger Rings(3pcs)
 - Cup Rings(2pcs)



NEW

DIA - ENDORING

- ITEM # 850 - 001
 - A convenient ring for Hand & Rotary files
 - Gauges from 1 to 30 mm
 - A simple & practical design for measuring files
 - A gold-plated ring



DIAFLEX™

- Flexible soft polypropylene tip syringe for root canal irrigation.
- Autoclavable and chemoclavable for multiple use up to 50 times.
- 1.3ml volume imprinted at intervals of 1mm.
- Provides maximum clinical effect in a minimum of time.


DiaFlex - III

- ITEM #1003-301
 - Package of 3 color coded syringes, 20 changeable tips and silicone oil


DiaFlex - V

- ITEM #1003-501
 - Package of 5 color coded syringes, 20 changeable tips and silicone oil

- ITEM #1001-315
 - Blue 20 disposable tips in a blister package



Tip Specification



INSTRUMENTS

PROCEDURE SET UP TRAYS & TUBES

STORAGE TUB(STD. SIZE-S)

- ITEM #719-101 White
 - ITEM #719-102 Gray
 - ITEM #719-103 Yellow
 - ITEM #719-104 Mauve
 - ITEM #719-105 Beige
 - ITEM #719-106 Blue
- Portable drawer for storage, organization, and transport of support materials used with procedure tray set-up
 - Dimensions : 323 x 57 x 280 m/m(W x H x D)



STORAGE TUB INSERT(STD. SIZE S)

- ITEM #719-201 White
 - ITEM #719-202 Gray
 - ITEM #719-203 Yellow
 - ITEM #719-204 Mauve
 - ITEM #719-205 Beige
 - ITEM #719-206 Blue
- Fits inside Storage Tub(Std. Size S) only, Sturdy 10 compartment design
 - Dimensions : 290 x 28 x 235 m/m(W x H x D)



PROCEDURE TUB

- ITEM #719-301 White
 - ITEM #719-302 Gray
 - ITEM #719-303 Yellow
 - ITEM #719-304 Mauve
 - ITEM #719-305 Beige
 - ITEM #719-306 Blue
- For storage & transportation of procedure materials
 - Dimensions : 310 x 70 x 275 m/m(W x H x D)



MINI FLAT TRAY(SIZE-F)

- ITEM #715-101 White
 - ITEM #715-103 Yellow
 - ITEM #715-106 Blue
- Sturdy, non-porous/stain-resistant polypropylene plastic trays
 - Available in flat
 - Dimensions : 243 x 22 x 168m/m(W x H x D)



B-LOCK FLAT TRAY

- ITEM #715-303 Yellow
 - ITEM #715-304 Mauve
 - ITEM #715-306 Blue
- Sturdy, non-porous/stain-resistant polypropylene plastic trays
 - Dimensions : 340 x 22 x 245m/m(W x H x D)



FLAT TUB SLIDE TRAY

- ITEM #715-601 White
- For storage & transportation of procedure materials
 - Dimensions : 244 x 28 x 133m/m(W x H x D)



DIVIDED TUB SLIDE TRAY

- ITEM #715-501 White
- Slides on top of dividers
 - Dimensions : 244 x 27 x 134m/m(W x H x D)



B-LOCK DIVIDED TRAY

- ITEM #715-203 Yellow
 - ITEM #715-204 Mauve
 - ITEM #715-206 Blue
- Sturdy, non-porous/stain-resistant polypropylene plastic trays
 - Dimensions : 338 x 22 x 245m/m(W x H x D)



INSTRUMENTS


INSTRUMENT CASSETTES

PLASTIC STERI CAGE(TYPE A)

- ITEM #714-103 Yellow
- ITEM #714-105 Beige
- ITEM #714-106 Blue



- Economical containers for processing procedure set-ups and individual instruments in ultrasonic cleaners and sterilizers
- Dimensions : 205 x 47 x 55m/m(W x H x D)




PLASTIC STERI CAGE(TYPE B)

- ITEM #714-203 Yellow
- ITEM #714-205 Beige
- ITEM #714-206 Blue



- Economical containers for processing procedure set-ups and individual instruments in ultrasonic cleaners and sterilizers
- Dimensions : 205 x 47 x 55m/m(W x H x D)



STERI COMPACT CASSETTE(SIZE A)

- ITEM #718-101 White
- ITEM #718-103 Yellow
- ITEM #718-106 Blue



- Steri-instrument cassette
- Holds 8 hand instruments or 16 when stacked
- Dimensions : 227 x 34 x 128m/m(W x H x D)



STERI COMPACT CASSETTE(SIZE B)

- ITEM #718-202 Gray
- ITEM #718-204 Mauve
- ITEM #718-206 Blue



- For processing oversized or bulky items
- Stackable
- Dimensions : 180 x 18 x 99m/m(W x H x D)



STERI SYSTEM CASSETTE(SIZE AA)

- ITEM #718-305 Beige



- Instrument holder included (capacity of 8 instruments)
- Stackable
- Dimensions : 185 x 35 x 140m/m(W x H x D)



STERI SYSTEM CASSETTE(SIZE BB)

- ITEM #718-405 Beige



- Fits standard-size ultrasonic cleaners
- Dimensions : 188 x 70 x 140m/m(W x H x D)




MINI INSTRU. MAT

- ITEM #702-102 Gray
- ITEM #702-106 Blue
- ITEM #702-107 Pink



- Holds up to 8 instruments
- Dimensions : 130 x 16 x 103m/m(W x H x D)




STERI INSTRU. MAT

- ITEM #702-202 Gray
- ITEM #702-206 Blue
- ITEM #702-207 Pink





- Holds up to 12 instruments
- Dimensions : 190 x 16 x 100m/m(W x H x D)



INSTRUMENTS



INSTRUMENT STORAGE BOX

DIVIDED TRAY(SIZE B)

- ITEM #715-403 Yellow
 - ITEM #715-404 Mauve
 - ITEM #715-406 Blue
- 
- 
- Sturdy, non-porous/stain-resistant polypropylene plastic trays
 - Available in divided style
 - Dimensions : 338 x 23 x 243m/m(W x H x D)






STAINLESS STEEL TRAY TYPE A

- ITEM #914-001
- 
- 
- A larger Tray for all purposes and storing.
 - Dimensions: 180 x 40 x 135m/m (W x H x D)






STAINLESS STEEL RACK TYPE A

- ITEM # 910-001
- 
- 
- Fits Most Sterilizing Systems.
 - Holds 5 Clip'in Instrument Holder Type A(See Next)
 - Dimensions: 190 x 95 x 97m/m (W x H x D)





STAINLESS STEEL RACK TYPE B

- ITEM # 910-002
- 
- 
- Fits Most Sterilizing Systems.
 - Holds 5 Clip'in Instrument Holder Type B(See Next)
 - Dimensions: 190 x 95 x 135m/m (W x H x D)





STAINLESS STEEL RACK TYPE C

- ITEM #910-003
- 
- 
- Holds 6 Sterilizing trays.
 - Compatible Stainless Steel DiaDent products include Clip'in Tray Type B, Instrument Cassette Type B, Stainless Steel Tray Type A.
 - Compatible Aluminum DiaDent Products include Container Type A(Sterilizing) & Instrument Tray Type A(Perforated & Non-Perforated).
 - Dimensions: 152 x 306 x 192m/m (W x H x D)



STAINLESS STEEL RACK TYPE D

- ITEM #910-004
- 
- 
- Holds 6 Sterilizing trays.
 - Compatible Stainless Steel DiaDent products include Instrument Cassette Type C, Stainless Steel Tray with Cover Type G.
 - Compatible Aluminum DiaDent Products include Instrument Tray Type B(Perforated & Non-Perforated).
 - Dimensions: 292 x 306 x 192m/m (W x H x D)



INSTRUMENTS

INSTRUMENT STORAGE BOX

COMPARTMENT BOX TYPE A

ITEM #822-001



- Holds 16 basins.
- Small Item Organizer.
- Dimensions: 230 x 25 x 110m/m (W x H x D)



COMPARTMENT BOX TYPE B

ITEM #822-002



- Holds 12 basins.
- Compartments vary in Size.
- Dimensions: 230 x 25 x 110m/m (W x H x D)



INSTRUMENT TRAY TYPE A (Perforated)

ITEM #826-001



- Holds 10 Instruments
- Ideal for Instrument Sterilization.
- Perforations Facilitate Sterilization & Drainage.
- Dimensions: 190 x 35 x 150m/m (W x H x D)



INSTRUMENT TRAY TYPE A (Non-Perforated)

ITEM #825-001



- Holds 10 Instruments
- Ideal for Instrument Storage.
- Dimensions: 190 x 35 x 150m/m (W x H x D)




INSTRUMENT TRAY TYPE B (Perforated)

ITEM #826-002



- Holds 16 Instruments
- Ideal for Instrument Sterilization.
- Perforations Facilitate Sterilization & Drainage.
- Dimensions: 280 x 35 x 180m/m (W x H x D)




INSTRUMENT TRAY TYPE B (Non-Perforated)

ITEM #825-002



- Holds 16 Instruments
- Ideal for Instrument Storage.
- Dimensions: 280 x 35 x 180m/m (W x H x D)




INSTRUMENT CASSETTE TYPE A

ITEM # 913-001



- Holds 14 Instruments.
- Has 4 silicone instrument holding units.
- Easy to open & close.
- Locking system.
- Dimensions: 215 x 25 x 105m/m (W x H x D)



INSTRUMENT CASSETTE TYPE B

ITEM #913-002



- Holds 20 Instruments.
- Has 4 silicone instrument holding units.
- Easy to open & close.
- Locking system.
- Fits Stainless Steel Rack Type C.
- Dimensions: 185 x 25 x 145m/m (W x H x D)



INSTRUMENTS


INSTRUMENT TRAYS

INSTRUMENT CASSETTE TYPE C

· ITEM # 913-003



- Holds 40 Instruments.
- Has 4 silicone instrument holding units.
- Easy to open & close.
- Locking system.
- Fits Stainless Steel Rack Type D.
- Dimensions: 280 x 25 x 190m/m (W x H x D)




PLASTIC CEMENT SPATULA

· ITEM #708-101 White



- High grade plastic that provide a sturdy handle and flexible tip for optimum mixing



CABINET TRAY



- CABINET TRAY #16
 - Hold 6 Instruments. Durable Plastic.
 - Dimensions: 195 x 20 x 95m/m (W x H x D)
- CABINET TRAY #16A
 - Hold 6 Instruments. Durable Plastic.
 - Dimensions: 195 x 20 x 95m/m (W x H x D)
- CABINET TRAY #17
 - 3 Compartment Plastic Tray.
 - Dimensions: 195 x 20 x 95m/m (W x H x D)
- CABINET TRAY #18
 - 3 Compartment Plastic Tray.
 - Long Lasting and Easy to Clean.
 - Commonly used in the operator.
 - Dimensions: 195 x 20 x 95m/m (W x H x D)
- CABINET TRAY #19
 - Multifunctional Plastic Tray.
 - Dimensions: 195 x 20 x 95m/m (W x H x D)
- CABINET TRAY #20
 - 6 Compartment Plastic Tray.
 - Dimensions: 195 x 20 x 95m/m (W x H x D)



#16



#16A



#17



#18



#19



#20

Color	#16	#16A	#17	#18	#19	#20
White	#713-101	#713-201	#713-301	#713-401	#713-501	#713-601
Yellow	#713-102	#713-202	#713-302	#713-402	#713-502	#713-602
Blue	#713-103	#713-203	#713-303	#713-403	#713-503	#713-603
Mauve	#713-104	#713-204	#713-304	#713-404	#713-504	#713-604
Beige	#713-105	#713-205	#713-305	#713-405	#713-505	#713-605
Gray	#713-106	#713-206	#713-306	#713-406	#713-506	#713-606

WAX KNIFE



Jujube wood handle knives

- Item #930-001 Size 1
- Item #930-002 Size 2
- Item #930-003 Size 3
- Item #930-004 Size 4
- Item #930-005 Size 5
- Item #930-006 Size 6
- Item #930-007 Size 7
- Item #930-101 Size Sharp
- Item #930-102 Size Round
- Item #930-103 Size M-type knives #134



M#134

round


sharp

#07

TRIPLE BITE IMPRESSION TRAY

- Beige Color.
- For accurate impressions and bite registrations.
- Thin non-tear mesh designed for extra retention of impression material.
- Rigid plastic frames and ergonomic Handles.
- Designed to minimize distortions.

- Anterior ITEM #720-001 : Pkg. of 35 units
- Posterior ITEM #720-002 : Pkg. of 48 units
- Quadrant ITEM # 720-003 : Pkg. of 40 units
- Sideless Posterior ITEM #720-004 : Pkg. of 48 units
- Assorted ITEM #720-005 : Pkg. of 40 units



MISCELLANEOUS ITEMS

DENTAL WAX

NEW

STICKY WAX


- ITEM # 940-010 Beige
- Made of superior resins and wax.
- Used for the temporary joint in the denture repair works.
- Easy to scrape off the denture after the plaster case is made.
- Diameter: 7m/m, Height: 80m/m
- 100g/box



NEW

INLAY BLOCK WAX

- ITEM # 940-030 Green
- Used for Inlay casting.
- Low congelation and contraction, and adequate cohesion and solidity levels make it possible to produce a more precise waxing pattern.
- Dimensions: 71 x 11 x 78m/m (W x H x D)
- 450g/box



NEW

WAX SEPARATOR

- ITEM # 940-050 Yellow
- Used for wax separation.
- Quick drying and not to stick to the denture surface, easily can be separated after some wax-ups.
- 2 bottles of 60ml and 1 bottle of 15ml



NEW

COLOR SPACER

- ITEM #940-060 SILVER
- ITEM #940-061 GOLD
- ITEM #940-062 RED

- One round of application generates space of 12 and 13µm.
- Before use, shake the bottle of Color Spacer.
- When the solution gets stiff, add a few drops of Diluent.
- 1 bottle of 60ml, 1 bottle of 15ml and 1 bottle(Diluent) of 60ml



SILVER



GOLD



RED

NEW

UTILITY WAX

- ITEM # 940-020 White
- ITEM # 940-021 Red
- Used to build a bank along the tray.
- Has the superior flexibility
- Diameter: 5m/m, Height: 315m/m
- 250g/box



NEW

MODELING WAX

- ITEM # 940-040 Pink
- Used for the denture works.
- Provided the sufficient stability for the base plate.
- The solidity remains equal both in the oral and room temperature.
- Dimensions: 170 x 1.4 x 100m/m (W x H x D)
- 500g/box



NEW

DIE HARDENER

- ITEM # 940-070 Blue
- Makes the surface of abutment stone harden to be consolidated and coated and prevent to be abrasive.
- 2 bottles of 60ml and 1 bottle of 15ml



MISCELLANEOUS ITEMS

NICKEL TITANIUM ARCH WIRES

NICKEL TITANIUM ARCH WIRES

- Last up to 5 times more than Stainless steel
- Larger working range than Stainless steel
- Less distortion and more shape memory than Stainless steel
- Easy to arrange in the bracket slots

SUPER ELASTIC ARCH WIRES

- Made of the highest quality Nickel-Titanium wire
- Highly polished ultra smooth finishing
- Allow a high degree of deflection
- Comes in Round(for class 1, mild cases) and Rectangular (for more stiffness and three dimensional control) Wires.

Round Type: Pkg. of 10

NT311-121	.012 Upper	SE-D-12-U
NT311-122	.012 Lower	SE-D-12-L
NT311-141	.014 Upper	SE-D-14-U
NT311-142	.014 Lower	SE-D-14-L
NT311-161	.016 Upper	SE-D-16-U
NT311-162	.016 Lower	SE-D-16-L
NT311-181	.018 Upper	SE-D-18-U
NT311-182	.018 Lower	SE-D-18-L
NT311-201	.020 Upper	SE-D-20-U
NT311-202	.020 Lower	SE-D-20-L

Rectangular Type: Pkg. of 10

NT321-661	.016 x .016 Upper	SE-R-1616-U
NT321-662	.016 x .016 Lower	SE-R-1616-L
NT321-621	.016 x .022 Upper	SE-R-1622-U
NT321-622	.016 x .022 Lower	SE-R-1622-L
NT321-721	.017 x .022 Upper	SE-R-1722-U
NT321-722	.017 x .022 Lower	SE-R-1722-L
NT321-751	.017 x .025 Upper	SE-R-1725-U
NT321-752	.017 x .025 Lower	SE-R-1725-L
NT321-821	.018 x .022 Upper	SE-R-1822-U
NT321-822	.018 x .022 Lower	SE-R-1822-L
NT321-851	.018 x .025 Upper	SE-R-1825-U
NT321-852	.018 x .025 Lower	SE-R-1825-L
NT321-951	.019 x .025 Upper	SE-R-1925-U
NT321-952	.019 x .025 Lower	SE-R-1925-L
NT321-151	.021 x .025 Upper	SE-R-2125-U
NT321-152	.021 x .025 Lower	SE-R-2125-L

REVERSE CURVE

Round Type: Pkg. of 10

NT313-121	.012 Upper	RC-D-12-U
NT313-122	.012 Lower	RC-D-12-L
NT313-141	.014 Upper	RC-D-14-U
NT313-142	.014 Lower	RC-D-14-L
NT313-161	.016 Upper	RC-D-16-U
NT313-162	.016 Lower	RC-D-16-L
NT313-181	.018 Upper	RC-D-18-U
NT313-182	.018 Lower	RC-D-18-L
NT313-201	.020 Upper	RC-D-20-U
NT313-202	.020 Lower	RC-D-20-L

Rectangular Type: Pkg. of 10

NT323-661	.016 x .016 Upper	RC-R-1616-U
NT323-662	.016 x .016 Lower	RC-R-1616-L
NT323-621	.016 x .022 Upper	RC-R-1622-U
NT323-622	.016 x .022 Lower	RC-R-1622-L
NT323-721	.017 x .022 Upper	RC-R-1722-U
NT323-722	.017 x .022 Lower	RC-R-1722-L
NT323-751	.017 x .025 Upper	RC-R-1725-U
NT323-752	.017 x .025 Lower	RC-R-1725-L
NT323-821	.018 x .022 Upper	RC-R-1822-U
NT323-822	.018 x .022 Lower	RC-R-1822-L
NT323-851	.018 x .025 Upper	RC-R-1825-U
NT323-852	.018 x .025 Lower	RC-R-1825-L
NT323-951	.019 x .025 Upper	RC-R-1925-U
NT323-952	.019 x .025 Lower	RC-R-1925-L
NT323-151	.021 x .025 Upper	RC-R-2125-U
NT323-152	.021 x .025 Lower	RC-R-2125-L

SUPER ELASTIC OPEN SPRING

Pkg. of 2

NT331-010	.010	SE-S-10-O
NT331-012	.012	SE-S-12-O

SUPER ELASTIC CLOSE SPRING

Pkg. of 10

NT332-010	.010	SE-S-10-C
NT332-012	.012	SE-S-12-C

ROUND WIRE

Pkg. of 5m

NT341-012	.012	RW-12-5
NT341-014	.014	RW-14-5
NT341-016	.016	RW-16-5
NT341-018	.018	RW-18-5
NT341-020	.020	RW-20-5

THERMAL ACTIVATED ARCH WIRES

- Achieve gentle torque force at the early stage of treatment
- Achieve optimum activation at 37.5°C (Mouth Temperature) and then the force level is increased
- Easy Manipulation

Round Type: Pkg. of 10

NT312-121	.012 Upper	SM-D-12-U
NT312-122	.012 Lower	SM-D-12-L
NT312-141	.014 Upper	SM-D-14-U
NT312-142	.014 Lower	SM-D-14-L
NT312-161	.016 Upper	SM-D-16-U
NT312-162	.016 Lower	SM-D-16-L
NT312-181	.018 Upper	SM-D-18-U
NT312-182	.018 Lower	SM-D-18-L
NT312-201	.020 Upper	SM-D-20-U
NT312-202	.020 Lower	SM-D-20-L

Rectangular Type: Pkg. of 10

NT322-661	.016 x .016 Upper	SM-R-1616-U
NT322-662	.016 x .016 Lower	SM-R-1616-L
NT322-621	.016 x .022 Upper	SM-R-1622-U
NT322-622	.016 x .022 Lower	SM-R-1622-L
NT322-721	.017 x .022 Upper	SM-R-1722-U
NT322-722	.017 x .022 Lower	SM-R-1722-L
NT322-751	.017 x .025 Upper	SM-R-1725-U
NT322-752	.017 x .025 Lower	SM-R-1725-L
NT322-821	.018 x .022 Upper	SM-R-1822-U
NT322-822	.018 x .022 Lower	SM-R-1822-L
NT322-851	.018 x .025 Upper	SM-R-1825-U
NT322-852	.018 x .025 Lower	SM-R-1825-L
NT322-951	.019 x .025 Upper	SM-R-1925-U
NT322-952	.019 x .025 Lower	SM-R-1925-L
NT322-151	.021 x .025 Upper	SM-R-2125-U
NT322-152	.021 x .025 Lower	SM-R-2125-L



MISCELLANEOUS ITEMS

GLASS ITEMS

GLASS DAPPEN DISHES



- ITEM #601-001 Clear
- ITEM #601-002 Amber
- ITEM #601-003 Blue
- ITEM #601-004 Green

- Standard Dappen Dish.
- Do not Autoclave.
- Easy to clean.
- Sold Individually or Box of 12pcs.



MEDICAMENT BOTTLES

- Locking lids for good sealing capabilities.
- Available in 4 colors for coding of chemicals.
- Easy to use and store.
- High quality glass, economically priced.
- 10ml volume
- Do not Autoclave.



ROUND TYPE

- ITEM #604-001 Clear
- ITEM #604-002 Amber
- ITEM #604-003 Blue
- ITEM #604-004 Green



SQUARE TYPE

- ITEM #605-001 Clear
- ITEM #605-002 Amber
- ITEM #605-003 Blue
- ITEM #605-004 Green



MIXING SLAB



- High quality, competitively priced Mixing Slabs available in 5 sizes.
- Do not Autoclave.



- ITEM #603-005 LT
- Dimensions : 150 x 15 x 80m/m (W x H x D)
- ITEM #603-004 LM
- Dimensions : 150 x 10 x 75m/m (W x H x D)
- ITEM #603-003 M
- Dimensions : 122 x 5 x 75m/m (W x H x D)
- ITEM #603-002 L
- Dimensions : 112 x 6 x 76m/m (W x H x D)
- ITEM #603-001 S
- Dimensions : 76 x 6 x 65m/m (W x H x D)

CHAIR SIDE TABLE

- ITEM #712-001 White



- Comes with attachments.
- Can be adapted to several types of chairs.
- Hole Diameter : available from 25 to 50m/m.
- Easy access to instruments.
- Dimensions : 390 x 30 x 235m/m (W x H x D)



GLASS ENDO ORGANIZER

- ITEM #602-000



- Shallow compartment walls for liquid sterilization.
- 6 inch diameter dish with 12 shallow wall compartments.
- Glass lid with depressions for chemicals.
- Do not Autoclave.



MEDICAMENT BOTTLE STAND

- ITEM #606-001



- Fits 5 Round or Square Medicament Bottles.
- Opaque high quality glass.
- Do not Autoclave.
- Color : Milky White
- Dimensions : 215 x 27 x 55m/m (W x H x D)



BIOHAZARD HANDY HOOPER

- ITEM #704-101 White



- A small plastic waste receptacle for biohazard items
- Can be wall mounted
- Dimensions : 170 x 160 x 100m/m (W x H x D)



TISSUE HOLDER

- ITEM #712-101 White



- Can be adapted to several types of chairs.
- Tissue box holder fits on unit chair pole.
- Hole Diameter : available from 25 to 50m/m.
- Dimensions : 245 x 54 x 125m/m (W x H x D)



MISCELLANEOUS ITEMS

RETRACTOR

X-RAY FILM HOLDER

- ITEM #750-001



- Used for Anterior and Posterior X-Ray.
- Angled for patient comfort in Anterior x-rays.



RETRACTOR (TYPE A)

- ITEM #740-001 Adult
- ITEM #740-002 Child



- Soft Flexible Plastic
- Pkg. of 2



SILICONE GRIPS

- Easy to make non-slip grip
- Available Color: Orange, Blue, Gray, Green & Brown
- Packed 5 different color Silicone Grips in poly bag




- ITEM #944-191 Size #1 (L:107mm, D:5mm)
- ITEM #944-291 Size #2 (L:88mm, D:5mm)
- ITEM #944-491 Size #4 (L:46mm, D:5mm)
- ITEM #944-791 Size #7 (L:42mm, D:5mm)



COLOR CODE CLIP

- ITEM #706-101 White
- ITEM #706-102 Gray
- ITEM #706-103 Yellow
- ITEM #706-104 Mauve
- ITEM #706-105 Beige
- ITEM #706-106 Blue



- 4pcs in a polybag
- Snap-in clips for trays, tubs and steri system cassettes
- Dimensions : 26 x 8 x 8m/m(W x H x D)

DAPPEN DISHES (PLASTIC)



- Color coded to help prevent cross contamination.

Color	50pcs	100pcs	250pcs
Clear	710-101	710-201	710-001
White	710-102	710-202	710-002
Yellow	710-103	710-203	710-003
Red	710-104	710-204	710-004
Blue	710-105	710-205	710-005
Black	710-106	710-206	710-006
Assorted	N.A	N.A	710-007



CHEEK RETRACTOR (TYPE B)

- ITEM #741-001 Adult
- ITEM #741-002 Child



- Soft Flexible Plastic
- Expands the mouth widely
- Pkg. of 2



ID-RING

- Silicone rubber
- Standard size(inner diameter 3.6mm)
- Pkg. of 120



- ITEM #942-101 White
- ITEM #942-102 Black
- ITEM #942-103 Brown
- ITEM #942-104 Green
- ITEM #942-105 Orange
- ITEM #942-106 Red
- ITEM #942-107 Gray
- ITEM #942-108 Yellow
- ITEM #942-109 Pink
- ITEM #942-110 Light Brown
- ITEM #942-191 Assorted



ENDO TRAINING BLOC

- ITEM #943-100



- Clear acrylic blocs.
- Ideal for training bloc for new techniques.
- Used for .04 tapered NiTi Rotary Files(curved type)
- Pkg. of 10

MISCELLANEOUS ITEMS

AMALGAM WELL

NEEDLE CAPPER STAND

- ITEM #705-101 White



- One-Handed uncapping
- Syringe is conveniently positioned
- Secured to table top
- Dimensions : 45 x 36 x 31m/m(W x H x D)



COTTON ROLL DISPENSER

- ITEM #770-001 White



- Holds approximately 60 Cotton Rolls.
- Features pull out drawer for easy cotton dispensing.
- Lid flips back for refilling.
- Non-slip bottom.
- Adjustable width for Cotton Roll in 20, 30 and 40m/m.
- Dimensions: 170 x 101 x 52m/m(W x H x D)



AMALGAM WELL



- Non slip Type with Rubber Base

- ITEM #819-001 Type A
 - Diameter: 30m/m, Height 20m/m
- ITEM #819-002 Type B
 - Diameter: 20m/m, Height 27m/m
- ITEM #819-003 Type C
 - Diameter: 20m/m, Height 30m/m
- ITEM #819-004 Type D
 - Diameter: 23m/m, Height 23m/m
- ITEM #819-005 Type E
 - Diameter: 31m/m, Height 23m/m



Type A



Type B



Type C



Type D



Type E

DISPOSABLE PLASTIC BOWL

- ITEM #707-101 White



- Rigid Mixing Bowls
- Pkg. of 5
- Diameter : 125m/m, Height : 85m/m



MIRACLE SUCTION

Innovative Multifunctional Mirrors with a Suction Function

- These mirrors features suction on the bottom and front sides to free up your handles and save your time.
- A major benefit fewer instruments in the oral cavity.
- Sure Grip autoclavable silicone handles.
- Fits into all suction unit on the dental chair.
- Available Adapters in 11m/m & 16m/m.
- ITEM #901-001 Regular Package 11m/m
- ITEM #901-003 Regular Package 16m/m
Contents: 5 autoclavable mirrors with detachable handle, 2 hoses, 2 adapters, O-rings and cleaning brush.
- ITEM #901-002 Small Package 11m/m
- ITEM #901-004 Small Package 16m/m
Contents: 2 autoclavable mirrors with detachable handle, 1 hose and 1 adapter.



Regular Package



Small Package

ISO 9001 ISO 13485



www.diadent europe.com

DiaDent Network

DiaDent Group International(Korea & International)

19-9 Songjung-Dong Heungduk-Gu Cheongju-si Choongchong Buk Do, Rep. of Korea

Tel)82-43-266-2315 Fax)82-43-262-8658

E-mail:diadent@chollian.net

DiaDent Group International Inc.(Canada & USA)

11-3871 North Fraser Way, Burnaby B.C., Canada V5J 5G6

Tel)1-604-451-8851 Fax)1-604-451-8865

E-mail:diadent@diadent.com

DiaDent Europe B.V (Europe)

Antennestraat 70, 1322 AS Almere, The Netherlands

Tel)31-36-549-8607 Fax)31-36-536-7317

E-mail:diadent@diadent europe.com