

MICROBRUSH® INTERNATIONAL
PRODUCT CATALOG



Microbrush®

Convenient. Consistent. Control.™



Micro-Applicators

Bristle Brush Applicators

Extended Reach Applicators

Adhesive Tip Applicators

Applicator Tips

Syringe Tips

Moisture Control

DESIGNED TO
ENHANCE EVERY
RESTORATIVE
PROCEDURE



Microbrush®

Convenient.
Consistent.
Control.™



PATIENT COMFORT



ABOUT MICROBRUSH®

MICROBRUSH® INTERNATIONAL

develops, manufactures and markets innovative disposable dental products focused on restorative dental procedures. Our products are designed to enhance everyday procedures and provide Convenient. Consistent. Control.™

MICROBRUSH® INTERNATIONAL

is committed to making the best dental products possible. All products are designed and developed in cooperation with leading dental material manufacturers and individual dentists to ensure outstanding product design and performance. For its high-quality products, the company has been awarded by REALITY and the Dental Advisor with the most prestigious recognitions earned only by the best dental products.

Microbrush manufactures its products in the USA, utilizing high-tech equipment managed under strict quality assurance programs, and is certified to ISO13485 for the production of medical devices. Microbrush strives to provide total user satisfaction and procedural success with every application.

COMPLETE LINE OF PRODUCTS



MICRO-APPLICATORS

NEW!

TRU™ BY MICROBRUSH®

Premium Micro-Applicator Series.....3

MICROBRUSH® PLUS

Dispenser Series.....4

MICROBRUSH®

Tube Series5

MICROBRUSH® X

Extended Reach Applicators.....6



MOISTURE CONTROL

NEW!

DRY TIPS®

Saliva Absorption Pads.....7



ADHESIVE TIP APPLICATORS

MICROSTIX™

Smart Hold & Original Hold.....8

BRISTLE BRUSH APPLICATORS

ULTRABRUSH®

Regular & Fine Size.....9



SYRINGE TIPS

FLOWTHRU®

Microbrush®-Style.....10

FLOWTHRU®

Standard-Style.....10



NEW!



DESIGNED FOR YOU.

Truly designed for you, TRU™ by Microbrush® is the first in Microbrush International's extensive line of disposable applicators to feature a double-bending tip, a longer 12-centimeter handle and an Ultrafine 0.5-millimeter tip size option.



- T TACTILE CONTROL**
12-centimeter handle with stay-put bending portion designed to provide maximum control
- R RELIABLE CONSISTENCY**
Advanced adhesive technology designed for greater consistency
- U ULTIMATE ACCESS**
Double-bending tip and new Ultrafine 0.5-millimeter head size designed for greater access

TRU™ BY MICROBRUSH® - PREMIUM MICRO-APPLICATOR SERIES



*Sold in individual tubes with 80 applicators per tube

REGULAR SIZE (2.0MM)

Blue..... 291159
Green..... 291160

FINE SIZE (1.5MM)

Yellow..... 291161
Pink..... 291162

SUPERFINE SIZE (1.0MM)

Gold..... 291163
Silver..... 291164

ULTRAFINE SIZE (0.5MM)

Teal..... 291165
Orange..... 291166

Microbrush BRAND *Plus*
DISPENSER SERIES



VIVID COLORS

**MULTIPLE AWARD-WINNING
APPLICATOR WITH
IMPROVED FEATURES**

- **STIFFER HEAD**
allows for scrubbing
- **STRONGER BENDING
PORTION**
for improved control

For precise applications in areas of limited access

- #1 micro-applicator worldwide specially packaged in disposable cartridges to fit a unique and convenient dispenser
- Multiple applicator colors are available to identify different materials during a procedure
- Three precise head sizes, Regular (2.0 mm), Fine (1.5 mm) and Superfine (1.0 mm), to match the size of the preparation
- Package contains 4 refill cartridges of 100 applicators each

DISPENSER SERIES - REFILL CARTRIDGE



Each cartridge contains
100 applicators. 4
cartridges per package



Kit contains 1 **sion™**
Dispenser plus 400
applicators

REGULAR SIZE (2.0MM)

Blue..... PR400BL
Peach..... PR400PE
Purple..... PR400PU
Green..... PR400GR
Assorted..... PR400

FINE SIZE (1.5MM)

Pink..... PF400PI
Yellow..... PF400YE
Assorted..... PF400

SUPERFINE SIZE (1.0MM)

White..... PSF400

DISPENSER KITS

(**sion™** Dispenser + Applicators)

REGULAR SIZE (2.0MM)

Assorted..... PR400-KIT

FINE SIZE (1.5MM)

Pink/Yellow..... PF400-KIT

SUPERFINE SIZE (1.0MM)

White..... PSF400-KIT

Microbrush[®] BRAND TUBE SERIES



VIVID COLORS

MULTIPLE AWARD-WINNING APPLICATOR WITH IMPROVED FEATURES


- **STIFFER HEAD**
allows for scrubbing
- **STRONGER BENDING PORTION**
for improved control

#1 micro-applicator worldwide packaged in practical tubes


- Multiple colors are available to identify multi-step procedures
- Microbrush[®] applicators can be used with materials such as cements, etchants, sealants, bonding agents, hemostatic solutions, and conditioners
- Hygienic Dispensing System - Shake out one applicator at a time with the flip top cap open



3 PRECISE HEAD SIZES:



Regular (2.0 mm)



Fine (1.5 mm)



Superfine (1.0 mm)



Each tube contains 100 applicators. 4 tubes per package

REGULAR SIZE (2.0MM)

Blue..... MRB400
Peach..... MRS400
Purple..... MRP400
Green..... MRG400
Assorted..... MRA400

FINE SIZE (1.5MM)

Pink..... MFP400
Yellow..... MFY400
Assorted..... MFA400

SUPERFINE SIZE (1.0MM)

White..... MSF400


Microbrush[®] BRAND EXTENDED REACH APPLICATORS



Extended, extra-thin applicator for confined spaces

- Designed with stiffer, non-linting fibers intended for scrubbing
- Designed for even placement of material to inner surfaces of tooth wall
- Tip delivers fractions of a drop
- Bending portion bendable up to 90°
- Use with **sion[™]** Dispenser

CAN BE USED FOR THE FOLLOWING:



Bonding
Posts



Inlay / Onlay
Cementation



Perio /
Subgingival

sion[™] Dispenser



1 sion[™] Dispenser
and 100 Applicators
Black.....MPX

- Antimicrobial properties built in to protect the dispenser
- Manufactured with sion[™] technology: Built-in silver ions penetrate the cell wall of bacteria and disrupt key cell functions, suppressing bacteria growth on the dispenser
- Specially designed to release one disposable applicator at a time
- Note: The dispenser's antimicrobial properties do not protect its users or others against bacteria, viruses, germs or other disease organisms
- Always use standard hygiene practices

REFILL CARTRIDGE



100 Applicators
Black.....PX100

NEW!

SALIVA ABSORPTION CAPACITY 30 TIMES ITS WEIGHT

Far outlasts cotton rolls and other absorbent pads,
and doesn't lose moisture when handled.



- Effective replacement for moisture control attributes of rubber dams.
- Complete coverage of parotid gland and entire buccal mucosa with polyethylene film for cutting protection.
- Facilitates work far back in oral cavity.
- No loose fibers or other residues left in mouth.
- Flexible; accommodates to cheek movements.
- Comfortably stays in place, yet is easily removed with water spray.

CAN BE USED FOR THE FOLLOWING:

- | | |
|-----------------------|-------------------------|
| • Impressions | • Composites |
| • Occlusion Grinding | • Cementation |
| • Orthodontics | • Fissure Sealing |
| • Paint-On Techniques | • Fluoride Applications |



LARGE SIZE
White..... 291542


SMALL SIZE
White..... 291543

*50 Dry Tips included in each package




Adhesive tip applicators designed to
handle delicate items


- Available in two holding strengths:
Smart Hold and Original Hold
- Smart Hold technology lets the practitioner determine the holding strength
- No adhesive residue left on surfaces
- Ideal for veneers, inlays / onlays, crowns and brackets



pick ...



place ...



twist .

Available in two holding strengths:

**ADVANCED ADHESIVE
TECHNOLOGY**



Smart Hold

Pressure-sensitive adhesive with Smart Hold technology facilitates the placement of delicate items with the advantage of easier release.



Original Hold

The strongest adhesive securely holds larger items with extra holding strength.




SMART HOLD
Green.....-STIX64G

ORIGINAL HOLD
Blue.....-STIX64B

64 adhesive applicators
per box


Storage: Shelf life 30 months. Storage Temperature: 12-28 C. / 54-82 F.
Keep *micro-Stix* in cool dry place.



Bristle brush applicators for spreading and smoothing materials

- Bends close to tip for improved control and reach
- Chemical-resistant fibers designed to hold up to the harshest dental materials
- Available in two sizes: UltraBrush® 2.0 (200 fibers) for better smoothing and spreading action and the UltraBrush® 1.0 (100 fibers) for more precise applications
- Reclosable cartridges

Better Coverage
Maximum Control



- Whitening
- Etchants



- Bonding Agents
- Sealants



- Desensitizers
- Adhesive / Resins

UltraBrush® Regular Size (200 fibers)



- DISPENSER KIT**
1 sion™ Dispenser and 100 Applicators
Green.....U2D
- REFILL CARTRIDGE**
2x100 Applicators
100 Yellow / 100 Green.....U2R200

UltraBrush® Fine Size (100 fibers)



- DISPENSER KIT**
1 sion™ Dispenser and 100 Applicators
Blue.....U1D
- REFILL CARTRIDGE**
2x100 Applicators
100 Orange / 100 Blue.....U1R200



- Bendable up to 90° without crimping
- Designed with emphasis on safety, ease of use, and patient comfort
 - Protective ring at base designed to reduce risk of cross-contamination
 - Secure thread lock
 - Plastic design reduces needle phobia
 - Available in Microbrush®-Style and Standard-Style
 - Package contains 100 dispensing tips




MICROBRUSH-STYLE		
Gauge	Color	Product No.
18	Pink	FM18-100
20	Yellow	FM20-100
22	Black	FM22-100
25	Blue	FM25-100




STANDARD-STYLE		
Gauge	Color	Product No.
18	Pink	FS18-100
20	Yellow	FS20-100
22	Black	FS22-100
25	Blue	FS25-100

SAVE MATERIAL

Less Waste



Flowthru™



Traditional Metal tip

More Waste

RECOMMENDED USAGE GUIDE *

Material	Gauge / Size		Color
Etchants	25	0.25MM	Blue
Sealants	22	0.41MM	Black
Flowable Composites	20	0.60MM	Yellow
Desensitizer	18	0.97MM	Pink

* Recommended gauge/size may vary depending on material viscosity



Microbrush

Convenient. Consistent. Control.™

1376 Cheyenne Avenue, Grafton WI 53024
microbrush.com

ISO 13485
Certified